

Product datasheet for RN200836

Dynlt1 (NM_031318) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dynlt1 (NM_031318) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Dynlt1
Synonyms:	AGS2; Tctex1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN200836 representing NM_031318 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGACTTCCAGGCCTCCGAGGAGACTGCATTTGTTGTGGATGAAGTGAGCAACATTGTAAGGAGG
CTATAGAAAGTGCCATCGGCGGTAAACGCCTACCAGCACAGCAAAGTCAACCAGTGGACCACTAATGTTGT
AGAACAGACTTTGAGCCAACTACCAAACCTGGGAAACCATTTAAATACATTGTGACCTGTGTGATCATG
CAGAAGAATGGTCTGGTTACACACCGCAAGTCTCTGCTTCTGGACAGCTCCACGGACGGGAGCTGCA
CAGTCCGATGGGAGAACAAGACCATGTACTGCATCGTCAGTGCCTTCGGACTGTCCAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_031318
Insert Size:	342 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_031318.1](#), [NP_112608.1](#)

RefSeq Size: 698 bp

RefSeq ORF: 342 bp

Locus ID: 83462

UniProt ID: [Q9Z336](#)

Cytogenetics: 1q11

Gene Summary: mediates vesicle transport from early to late endosomes in cooperation with dynein [RGD, Feb 2006]