

Product datasheet for **RN200935**

Cutc (NM_001108525) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cutc (NM_001108525) Rat Untagged Clone
Tag: Tag Free
Symbol: Cutc
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN200935 representing NM_001108525
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGAAAGGGATTTCTTGGCAGCGGAAACAGGCGTGGATATCGTCGGGGAAGGCGGGAGCAGGAAATG
GATTTCTTATGGAAGTGTGCGTAGATTCGGTGGAGTCGGCTGTGAATGCAGAAAGAGGAGGTGCTGACCG
GATTGAATTATGTTCTGGTTTACTGGAGGGAGGAACACACCCAGCATGGGTGTTCTTCAAGTAGTAAAG
CAGAGTGTTCCAGATCCCAGTTTTGTGATGATTCGACCACGTGGAGGTGATTTTCTATATTCAGATCGTG
AAGTTGAGGTGATGAAGGCTGACATTCGTTGCCAAGCTTTATGGTCTGATGGTTTGGTATTTGGGGC
ATTGACTGAAGATGGACACATTGACCAAGAGCTATGCTTGTCTCTTGTAGCTCTCTGCCGTCCTTACCA
GTCACCTTCCACAGAGATCCAATGGCAGCTCTGGAGACCCTCCTGACCTTGGGGTTTGAACGCGTTCTGA
CCAGTGGGTGTGACAGTTCAGCACTGGAAGGACTGCCTCTCATAAAGCAGCTCACAGACCAGGCAAAAGG
CAGGATTGTGGTAATGCCAGGAGGTGGGATAACAGTTAAAAATCTGCAAAGGATCCTTGAGGGTTCTGGA
GCTACAGAATCCACTGTTCTGCCCGATCTTCTAGAGACTCAGGAATGAAGTTTCGAAATTCCTCTGTTG
CAATGGGAGCTTCACTTGCTCATTGAGAAATTTCTCTAAAGGTAACAGATGTGACAAAAGTAAGGACTTT
GAATGCTATTGCAAAGGATGTCCTAGTT**TAG**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_001108525
Insert Size: 801 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108525.1</u> , <u>NP_001101995.1</u>
RefSeq Size:	1302 bp
RefSeq ORF:	801 bp
Locus ID:	361760
Cytogenetics:	1q54