

## Product datasheet for **RN200939**

### Dgat2l6 (NM\_001109370) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dgat2l6 (NM\_001109370) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dgat2l6  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN200939 representing NM\_001109370  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCTTCTTCTCCCGTCTTGATCTGCATGAAGGCCTCCGAACACTCTCCGTATTACAGTGGATTCCAG  
TTTATGTAATTTAGGAACCTCTTCCATCTTGGGGATACCTATTTCTGCTATTCAGTACCTTTTGCC  
CTTGTCTGTGCTCTTCTAGCCTGGGTGCCTATGATTGGAACACTCACAGTCAAGATGGTAGGCCGTTCA  
GCTTGGGTACGAACTGGACCTCTGGAAGTATTTCAAAGCTACTTCCAGTAAAGCTGGTGAAGACCC  
ATGACCTTTCTCCAAGCACAACTACATCATTCTCAGCCACCCCATGGCATTCTCTGTTATGGTGCCT  
CATCAACTTTGCTACTGAAGCCACTGGCTTTTCCAGGGTTTTTCCCTTCTATCACTCCCTTTTAGCAACC  
TTAGAAGGAATCTTCTGGATCCCATTTGTTCGAGATTATGTGATGTCCATGGGTATTTGTCCAGTGAGTG  
AGTTGGCCCTGAAACACCAGCTGACCCAAAGAGGCTCAGGCAATGCTGTGATCATTGTGCCAGGTGGAGC  
TTCCGAAGCCCTCTTATGCTGCCAGGAGTCTCCACGGTCTACCTCAAAAAACGAAAAGTTTTGTGAAG  
CTAGCACTGAAGACAGGAGCATACTTGTCCCTTCTTACTCCTTTGGTGAGAATGAGGCTTACAACCAAG  
AGACCTTTCTGAAGGCACATGGCTAAGGTTCTTCAAAAAAAGTTCCAGAAAATGGCAAAAAGAACTCT  
AGGCATAAATTTCTGTACTTTTCATGGCAGGGGCTCACTAGGGGATCTGGGGCTTTTGCCTTCAAC  
CATCCCATTACCACTGTTGTTGGAGAGCCCTTCCAGTCCCAAGATCCTTGATCCTGACAAGGAGACCG  
TGGCAAAGTACTTTGAACTTTACATCAGTGCCTACGCAAGCTGTTTGACCAACACAAAAGTAAATATGA  
CTTCTCCAAGACTCAGGAGTTGACAATTATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001109370



[View online »](#)

<b>Insert Size:</b>	1014 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001109370.1</a></u> , <u><a href="#">NP_001102840.1</a></u>
<b>RefSeq Size:</b>	1386 bp
<b>RefSeq ORF:</b>	1014 bp
<b>Locus ID:</b>	678749
<b>Cytogenetics:</b>	Xq22