

## Product datasheet for **RN200973**

### Fundc2 (NM\_001135757) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fundc2 (NM_001135757) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fundc2
Synonyms:	RGD1560916
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN200973 representing NM_001135757 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGCCGCCAACGGTCAAGGGGACCCGACAGGAAGCGCGAGGCCCTGGACCTCGCGGAATTCACGCGGA  
GGCAACCATGGTGGCGTAAGCTGTTGGGCCAGAGTCTGGGCCCTCGGCTGAGAAGTACAGCGCGGCGAC  
GCAGCTGGCCGTCGGGGGTGTCACTGGATGGTGCACGGGGCTCGTGTCCGGAAGGTCGGAAGCTGGCA  
GCCACGGCCCTAGGAGGAGGGTTCTTCTCCTGCAGCTGGCAAGCCACACGGGGTACATCAGAGTCGACT  
GGCAGAGACTGGAAAAAGACGTAAGAAAGCCAAGGCGCAGCTGAAGATCCATAGGGGTCGGCAAATGCC  
CACTGAGGTGGGAGCAAAGCCGAGGAGGCCCTGTCGTTCTGTAAGAAGAATGTTCCAGTAGCCGAGGG  
TTTCTTGAGGCTTCTGCTCGGCCTGGCGTCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001135757
Insert Size:	456 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001135757.1, NP_001129229.1</u>
<b>RefSeq Size:</b>	1397 bp
<b>RefSeq ORF:</b>	456 bp
<b>Locus ID:</b>	361288
<b>Cytogenetics:</b>	18p13