

## Product datasheet for **RN201009**

### Cyp39a1 (NM\_001106893) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cyp39a1 (NM\_001106893) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cyp39a1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN201009 representing NM\_001106893  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGCATTATGGAGTTATTCTACCAATAGCAATAGCCGTCCTAGGAAGCTGCGTTCTTTTCTATTCT  
CTCGGTGGAAGAATTTGCGTGGACCCCGTGCATCCAGGGTTGGATCCCTGGATCGGAGCTGGCTTTGA  
GTTTGGGAAAGCCCTCTTGAATTTATAGAAAAAGCAAGAATCAAGTATGGACCAGTGTACTGTCTTT  
GCTGTGGGAAAACGAATGACCTTTGTTACTGAAGAAGAAGGAATCAATGTCTTCTAAAGTCTAAGCATG  
TAGATTTTGAATTAGCAGTACAACGCCCTCTCTACCACACAGCATGGATCCCAAAGAATATTTTCTCGC  
ACTGCATGAAAACTGTATGTCTTGATGAAAGGGAAGTGGAACTTTCAACACACATCACTTTACTGGG  
CAGCTGACTGAAGAATTCATGATCAACTGGAAGGTTTAGGCACTCATGGGACTATGGACCTGAATGACT  
TTGTACGGTATCTACTTTATCCAGCCACACTCAATACTCTCTTTATGAAAGGGCTGTTTCTCACAGACAA  
GAGGAAGATCAAGGAGTCTACCAGCATTTTAAACCTATGATGAGGGCTTTGAGTATGGCTCCCAGCTG  
CCCGAATGGCTTCTGAGAACTGGTCAAATCCAAAAGATGGCTCTTGGCACTGTTTGAAGAAAATATTG  
GAGACATAAAAACTCATGGATCTGCAGGTCATTCCGAAACCTTTTGCAGGCCGTCTGGGAATGGTGGAA  
GACAGAGACACGTCTACACAGTCCAAATACGGACTGGTATGCTTTGGGCTTCTCTGGCCAATGCTGCT  
CCTATTGCGTTCTGGACACTTGCATATATCCTGTCTCATCCCGATCTGCACAGGACAATTGTAGAAAGCA  
TATCTTCTGTGTTTGGCACTGCAGGTAAGACAAGATTCAAGTGTCTGAAAATGACTTGAAGAAGCTTCT  
TCTAATTAAGTGGTGTATTCTGGAATCCATCCGCTTGAGAGCCCTGGGGTCATTACTAGGAAGGTAGTG  
AAGCCAGTTAAATTTCTGAATCACACAGTTCCTTCTGGAGATCTGCTGATGCTGTCTCCATTCTGGCTAC  
ACAGAAATCCAAAATATTTCCCTGAACCCAGAGTCTTTCAAACCTGAACGCTGGAAAGAGGCAAATTTAGA  
TAAGTACATTTTCTGGACTACTCATGGCTTTTCGGAGGCGGGAAGTTTCAGTGTCTCTGGAAGGTGGTTT  
GCTCTCTTAGAGATTCAACTGTGATTATCTTAGTGTCTTATAAGTATGAGTGCAGTCTTCTGGATCCAT  
TACCAAAGCAGGCAAAGTTCTCTTCATTTGGTGGGTGTCCCCAGCCAGCAGGAAATGCCGAAT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001106893
<b>Insert Size:</b>	1398 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001106893.1</a>, <a href="#">NP_001100363.1</a></u>
<b>RefSeq Size:</b>	3195 bp
<b>RefSeq ORF:</b>	1398 bp
<b>Locus ID:</b>	301264
<b>Cytogenetics:</b>	9q13