

Product datasheet for **RN201022**

Cbwd1 (NM_133535) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cbwd1 (NM_133535) Rat Untagged Clone
Tag: Tag Free
Symbol: Cbwd1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN201022 representing NM_133535
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTACCTGCTATGAAAACCGTGGAGGCGGAAGAAGAGTACGCGGAGGATTGTCCCGAATTGGTCCCCA
TTGAGACCAAGCACCAAGAAAAAGAGAAAATCTTGACTTCATCATTAAAGATTCCAGTCACAATTGTCAC
CGGGTATTTAGGTGCTGGAAAGACAACACTTCTGAATTACATTCTGACAGAACAACATAATAGGAAAAA
GCAGTGATCTTAAATGAATTTGGAGAAGGGAGTGCAGTAGAGAAATCCTTAGCTGTCAGTCAAGGAGGAG
AACTCTATGAAGAATGGCTGGAACCTAGAAACGGTTGCCTCTGCTGTTCCGGTCAAGGACAATGGCCTTAA
AGCTATAGAGAATCTGATGCAGAAGAAGGGCAAGTTTATTACATACTCTTAGAGACCACTGGGCTAGCA
GATCCTGGTGCAGTGGCTTCTATGTTTTGGGTTGATGCTGAGTTAGGAAGTGATATTTATCTTGATGGCA
TCATAACTGTTGTAGATTCAAAATATGGATTAAAAATTTAACTGAAGAAAAACCTGATGGCCTTGTCAA
TGAAGCTACTAGGCAAGTTGCTTTGGCAGACATGATTCTCATCAATAAAACAGACTTGGTTTCAGAAGAG
GAGCTAAACAAATGAGAACAACATCAGATCAATAAATGGACTGGGAAAAGTTTTAGAACTCAAAGAT
CAAGAACTCATCTCTCAAATATATTAGATCTTCATGCCTATGACACTCTCTCTGGTATAAGTTTGCAGAA
AAAACCTCCAGCATGTGTCAACGGCACCCACCTTGATCAGAGCATTGTTACAGTCACATTTGACGTACCA
GGAAGTGCAGAGGAAGAGAGTCTAAATGTATTTATCCAGAATCTCCTGTGGGAAAAGAATGTGAAGAACA
AGGACGGCCGCTGCATGGAGGTCATCCGGCTGAAGGGCCTGGTGTCTATCAAAGACAAACCACAGCAAAT
GATTGTGCAAGGCATCCATGAGCTGTACGAGCTGGAAGAGAGCCGAGTGAAGTGAAGGACGATGCCGAG
AGAGCCTGCCAGCTGGTCTTTATTGGTAAAAATTTAGATAAGGATATTCTTCAACAACTCTTCATCACTG
CTGTGGCAGAGACAAAAGAACAACGACAGCACCAGGCCAAGATGGAGCCTGTCCGCCGACT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_133535



Insert Size:	1185 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_133535.2</u> , <u>NP_598219.2</u>
RefSeq Size:	1405 bp
RefSeq ORF:	1185 bp
Locus ID:	171057
UniProt ID:	<u>Q99MB4</u>
Cytogenetics:	1q51