

Product datasheet for RN201034

Ccr4 (NM_133532) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccr4 (NM_133532) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ccr4
Synonyms:	Cmkbr4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201034 representing NM_133532 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATGCCACAGATGTCACAGACACCACCCAGGATGAAGCCGCGTACAATAGTTACTACTTCGACGAGA
GCTTGCCAAAGCCTTGCACCAAGGAAGGCATCAAGGCGTTCGGCGAGGTCTTCTGCCTCCTCTACTC
CTTGGTCTTTCTGTTGGGCTGTTTGGGAATTCTGTTGTGGTTCTGGTGTGTTCAAGTACAAGAGGCTC
AAGTCCATGACTGACGTGTACCTGCTGAACCTAGCCATCTCGGATTTACTGTTTGTGCTGTCTCTCCGT
TCTGGGGTACTATGCTGCTGACCAGTGGGTTTTGGTCTAGGTCTGTGAAGATTATCTCATGGATGTA
CCTGGTGGGCTTCTACAGCGGGATCTTCTTCATCATGCTCATGAGTATAGACAGATACCTGGCCATCGTG
CACGCAGTGTTCCTTGAGAGCGAGGACTCTGACCTATGGGGTCATCACCAGCCTCATCAGTGGTCTG
TGGCTGTGTTGCCTCCCTACCAGGCCTCTTGTTCAGCACCTGCGACACAGAGAACAACCACATACTG
CAAACCAGTACTCGGTCAACTCGACCAGTGGAAAGTCCTCAGCTCCCTGGAGATCAATGTCTTGGGG
CTGGTCATCCCCCTGGGATCATGTTGTTCTGCTATTCCATGATCATCAGGACCTGAGACACTGCAAGA
ATGAGAAGAAGAACAGAGCAGTTCGGATGATCTTCGCCGTGGTGGTCTCTTCTCGGCTTCTGGACGCC
GTACAACGTGGTGTCTTCTCGGAGACCTTGGTGGAGCTCGAAGTCCCTCAGGACTGCACCTTGGAGAGG
TACCTAGACTATGCCATCCAGGCCACAGAAACCCTGGCCTTCATCACTGCTGCCTTAACCCTGTCAATT
ACTTCTTCTTGGGAGAAATTCGCAAGTACATCGCCCAACTCTCAGAACGTGCCGGGTCTCTCTGT
GCTCTGCCAACACTGTGACTTCTTTCAGATCTCTCGGCCGACACTTCCAGCTCCTTACACTCAGTCC
ACTGTGGATCATGACTTCCGTGACGCTTTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2339_f06.zip



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_133532
Insert Size:	1083 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_133532.2</u> , <u>NP_598216.2</u>
RefSeq Size:	1083 bp
RefSeq ORF:	1083 bp
Locus ID:	171054
Cytogenetics:	8q32
Gene Summary:	mouse homolog is a chemokine receptor that binds both macrophage inflammatory protein-1 alpha and RANTES [RGD, Feb 2006]