

Product datasheet for **RN201149**

Crx (NM_021855) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Crx (NM_021855) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Crx
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Fully Sequenced ORF: >RN201149 representing NM_021855
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATGGCGTATATGAACCCGGGACCTCACTATTCGGTCAATGCCTTGGCCCTGAGTGGCCCCAGTGTGG
 ACCTGATGCACCAGGCTGTCCATACTCAAGTGCCCTAGGAAGCAGAGGCGGGAGCGGACCACGTTAC
 CAGGAGCCAGCTGGAGGAGCTGGAGGCCCTGTTTGCCAAGACCCAGTACCCGGATGTGTATGCACGTGAG
 GAGGTTGCTCTCAAGATCAACCTGCCTGAGTCCAGGGTCCAGGTCTGGTTCAAGAATCGCAGGGCGAAAT
 GCAGACGACAGAGCAGCAGCAGAAACAGCAGCAGCAGCCCCAGGGGTACAAGCCAAGGCTCGTCCTGC
 GAAGAGGAAGGCAGGCACATCCCCGAGACCCTCTACAGATGTTTGTACAGATCCTTTGGGCATCTCAGAT
 TCCTACAGCCCATCTCTGCCTGGTCCCTCAGGTTACCTACCACAGCAGTGGCCACCGTATCCATCTGGA
 GTCCGGCCTCAGAGTCTCCTTTGCCTGAGGCCAGAGAGCTGGCCTAGTGGCTTCAGGGTCTCTCTCAC
 CTCAGCCCCTTATGCCATGACCTATGCCCTGCTTCTGCTTTCTGTTCTTACCTTCAGCTTATGCATCT
 CCAAGCTCCTATTTAGTGGTTGGATCCATACCTTTCTCCATGGTGCCCGAGCTAGGGGGTCCGGCTC
 TCAGCCCCCTCTCGGGCCCTCTGTGGGGCCATCCTTGCCCGAGTCCCCACCTCCTTGTGAGCCAGAG
 CTATAGCACCTACAGCCCCGTGGACAGCCTGGAATTCAGGACCCTACAGGCACCTGGAATTCACCTAC
 AATCCCATGGACCCTCTGGACTACAAGATCAGAGTGCCTGGAAGTTTCAGATCTT**GTAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
 ACCN: NM_021855
 Insert Size: 900 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_021855.1</u> , <u>NP_068627.1</u>
RefSeq Size:	900 bp
RefSeq ORF:	900 bp
Locus ID:	60446
Cytogenetics:	1q21
Gene Summary:	an OTX-like homeobox protein; may play an important role in the retina [RGD, Feb 2006]