

Product datasheet for RN201158

PrI2b1 (NM_138861) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PrI2b1 (NM_138861) Rat Untagged Clone
Tag:	Tag Free
Symbol:	PrI2b1
Synonyms:	Ghd11; Prlpk
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201158 representing NM_138861 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGCTACTGTATTTGCCTCAAATATTTTCTCAAGGGCATCATCGCTCCTGTTCTGGTGCCTACTTGC
TCTTTTGGGAGAATGTAGCATCTATATCCACCTGTGCAGAGAGGGATGCCACCATTGAGCATTCTTAGA
GAACTTCTTACATTGACAACCTTTATGTCTCATGTGATGATTGAACTGCAAACTCTTCACTGAA
TTTAATAATCAGTATGCCAGGGTAAGAGATACAATGACAGGATCCCTGGAACGTGTCACTGCTTTT
TTGATACTCCAGTAAACAAGGAGCAATCTCTAGGAAGTGATCCGAAAACACTACTGAAATTTGGTACGCAG
TTTATTGAATTCCTGGACCAATGCTCTAAATCATCTTGTGAATGAAATATCTGCAATGCAAGGAGACCT
TCTTTTCTTTCTCCAAAGCCAGAGAGATTGAGCAAAATTTGATGAACTCATCGACGGGTGTAACAA
TTCTCAGCATGATTGGAGAGAGAGATAATGACACCCACCTTGCTTGGTCCGGGCTGTATCCTTGCAGTC
AAGCAATGAAGATGTTTCGCTGCTTTTCTTTTATACCTGATTGCTGCTGCTCCGAGACTCTCGAAAA
GTAATACTTATCTTGAGGTTATCAAGTACCAATATTCAATCAAACAACCTGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_138861
Insert Size:	687 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_138861.1</u> , <u>NP_620216.1</u>
RefSeq Size:	865 bp
RefSeq ORF:	687 bp
Locus ID:	192224
UniProt ID:	<u>Q9JKL9</u>
Cytogenetics:	17p12
Gene Summary:	a hormone that is involved in the control of maternal-fetal adaptations to pregnancy [RGD, Feb 2006]