

## Product datasheet for **RN201162**

### Gpr87 (NM\_001107677) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gpr87 (NM\_001107677) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gpr87  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN201162 representing NM\_001107677  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGC**C

ATGAATGTCTCTCCTTTTCAGGTAATGAATTATACAGCCAAGCAAGTCACACTGCAAACAGCAGCATCG  
AGGGACATGGGAAGAAGTCCACCCATCACAACGAGTTTGACACCATCATCCTGCCGGTGCTTTACCTAGT  
TATATTTGTGGCAAGCATCTGCTGAACGGTCTGGCTGTGTGGATCTTCCACATTCGGAATAAAACC  
AGCTTCATATTTTATCTCAAAAATATAGTAGTAGCTGACCTCATAATGACCCTGACATTCGGTTCCGAA  
TAGTCCGTGACGCAGGATTTCGGACCTTGGTACTTTGAGTTTATCCTCTGCAGATACACCTCGTTTTGT  
CTATGCAAACATGTATACATCCATTGTGTTTCTGGGCTGATCAGTGTTGATCGTATCTGAAGGTGGTA  
AAGCCTTTTGGTACTCTCGCATGTATAGCATAACCTTCAACAAAGTCTTATCGATTTGTGTCTGGTGA  
TCATGGCTGTTCTGCTTGCACAAATCATACTGACTAACGGTCAACCAACTAAGGAAAACATTCATGA  
CTGCATGAAGCTTAAAAGTCCGTTGGGAGCCAAGTGGGATATGGCTGTACCTACGTGGACAGCTGTCTG  
TTTGTGCTGTGCTGGTATTCTGATTGGATGCTACATAGCCATCTCCAGATACATCCACAAATCCAGCA  
GGCAATTCATAAGTCAATCGAGCCGGAAGCGGAAGCATAACCAGAGCATCCGCGTGGTCTGGCTGTGTT  
TTTTACCTGCTTCTCCGTATCACTTGTGAGAATCCCCTTCACTTTCACTTAGACAGGCTCCTA  
GATGACTCAGCACAAAAATCTCCACTACTGCAAGGAAATGACACTTTTCTGTCTGCGTGAACGCTGT  
GCCTGGATCCGATAATTTATTTTTTTCATGTGTAAGTCATTTTCAAGAAGTTATTCAAGAAATCAACAT  
AAGAACCAGGAGTAAAAGCATCAGGTCGCTGCAAAGCGTCCGAAGATCAGAAGTGCATCTACTATGAC  
TACACCGATGT**GTAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001107677  
**Insert Size:** 1065 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001107677.1</a></u> , <u><a href="#">NP_001101147.1</a></u>
<b>RefSeq Size:</b>	1161 bp
<b>RefSeq ORF:</b>	1065 bp
<b>Locus ID:</b>	310443
<b>Cytogenetics:</b>	2q26