

Product datasheet for **RN201180**

Adcyap1 (NM_016989) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adcyap1 (NM_016989) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Adcyap1
Synonyms:	Pacap
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201180 representing NM_016989 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCATGTGTAGCGGAGCAAGGTTGGCCCTGTTGGTCTACGGGATAATAATGCATAACAGCGTCTCCT
GTTACCTGCCGCCGACTCAGCTTCCCTGGGATCAGACCAGAAGAAGAGGCTTACGATCAGGACGGAAA
CCCGCTGCAAGACTTCTACGACTGGGACCCTCCGGCGCAGGGAGCCCCGCCTCCGCGTGCCTGACGCC
TACGCCCTTACTACCCAGCCGACAGGAGAGATGTCGCCACGAAATCCTAACGAAGCCTACCGCAAAG
TCTTGGACCAGCTGTCCGCCAGGAAGTACCTGCAGTCCATGGTGGCCAGGGGCATGGGCGAGAACCTCGC
CGCCGCCGCGGTGGACGACCGGGCACCCCTTACCAAACGCCACTCGGACGGCATCTTACAGACAGCTAT
AGCCGCTACCGAAAACAAATGGCTGTCAAGAAATACTTGGCGGCCGTGCTAGGGAAAAGGTATAAACAGA
GGGTTAAAAACAAAGGACGCCGAATAGCGTACTTGT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms:	https://cdn.origene.com/chromatograms/ja1923_c11.zip
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_016989
Insert Size:	528 bp



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016989.2 , NP_058685.1
RefSeq Size:	1262 bp
RefSeq ORF:	528 bp
Locus ID:	24166
UniProt ID:	P13589
Cytogenetics:	9q38
Gene Summary:	regulatory protein that stimulates cell proliferation, induces differentiation, modulates neurotransmitter release, vasodilation, bronchodilation, and hormone release, and suppresses inflammatory responses [RGD, Feb 2006]