

## Product datasheet for RN201227

### Kcnj13 (NM\_053608) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kcnj13 (NM\_053608) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Kcnj13  
**Synonyms:** Kir7.1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN201227 representing NM\_053608  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGACAGCAGGAATTGTAAGTTAATGCTCCTCTCCTAAGTCAAAGATACCGGAGGATGGTCACAAAAG  
ATGGCCACAGCACACTTCAAATGGATGGTCTCAAAGAGGTCTTGTGTATCTTCGGGATGCATGGGGAAT  
CCTAATGGACATGCGCTGGCGTGGATGATGCTGGTCTTTCTGCTTCTTCGTTGCCACTGGCTTGT  
TTTGACAGTCTCTGGTATGCTGTAGCTGAGATGAATGGTATCTGGAAATAGACCACGACGTCCCGCCTG  
AAAACCACACTATCTGTGTGAAGCACATACCAGTTTCACAGCTGCGTTCTCCTTCTCTCTGGAGACACA  
ACTTACAATCGGCTATGGCACCATGTTCCCCAGTGGGGACTGTCCAAGTGAATCGCCTTACTTGCATA  
CAAATGCTCCTAGGCCTCATGCTAGAGGCCCTTATCACAGGTGCCTTTGTGGCGAAGATTGCACGGCCAA  
AAAATAGAGCTTTCTCAATCCGATTCACAGACTTAGCAGTTGTAGCTCACAAAGATGGCAAACCTAATCT  
TATCTTCCAAGTGGCCAAACCCGGCCAGCCCTCTGACCAGCGTTCGGGTCTCTGCTGTACTCTATCAG  
GAAAGAGAAAATGGTGAAGTCTACCAGACCAGTGGACTTCCACCTTGATGGCATCAGTCTGAGGAAT  
GTCCATTCTCATCTTCCCGCTAACCTATTATCATACTATCACACCATCAAGTCCCTAGCTACCCTGCT  
CCAGCATGAGACTCCTTACACTTTGAATTAGTTGTTCCTCTCAGCAATGCAAGAAGGCACTGGGGAA  
ATCTGCCAAAGGAGGACATCCTACCTACCTCTGAGATCATGTTACATCACCGATTTGCAGCTCTAATGA  
CTCGAGGTTCCAAAGGTGAGTATCAAGTCAAGATGGAGAATTTTGACAAGACTGTTCTGAACATCCAAC  
TCCTGTGGTCTCTAAGAGTCCACATAGGACTGACCTAGATATTATATCAATGGACAAAGCATTGACAAT  
TTTCAGATTGCTGAAACAGGGCTGACAGAA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_053608



<b>Insert Size:</b>	1083 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_053608.2</a></u> , <u><a href="#">NP_446060.1</a></u>
<b>RefSeq Size:</b>	1268 bp
<b>RefSeq ORF:</b>	1083 bp
<b>Locus ID:</b>	94341
<b>UniProt ID:</b>	<u><a href="#">O70617</a></u>
<b>Cytogenetics:</b>	9q35
<b>Gene Summary:</b>	inwardly rectifying K <sup>+</sup> channel; may be involved in recycling K <sup>+</sup> for the Na <sup>+</sup> ,K <sup>(+)</sup> -ATPase [RGD, Feb 2006]