

Product datasheet for RN201260

Uts2b (NM_198133) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Uts2b (NM_198133) Rat Untagged Clone
Tag: Tag Free
Symbol: Uts2b
Synonyms: U2B; Urp; Uts2d
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN201260 representing NM_198133
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGTTCTTCTCAACTTCTCTGCTTTGGACTCCTGGCTTTGTTGTCTGTGACGACCCTTTACATT
CTGTACGTGGACGGCCGATCTGAGCTCAGGACATGAGTTATTTCCAGCTGAAGAACACAGACTCAGGA
GAAACTACCGCTGGTCTGCTGATCCGAAACCCTGGTTTCCAGAGCCCGCCCATGCTGGTGTAGATTTG
CCAGCAAAGTGAAGAAGCTTAGACAGCTAAAGAAGCTGAGAGAATGGTTTATGGAGCAAAGAGTGCTG
AGCCGTCCAATGCTTTAGACAAGCTATCCTCTCCACCCTATTAACGAGCTTGCTTTGGAAACTG
TGTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_198133

Insert Size: 357 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198133.2](#), [NP_937766.1](#)

RefSeq Size: 747 bp

RefSeq ORF: 357 bp

Locus ID: 378939

UniProt ID: [Q76512](#)

Cytogenetics: 11q22

Gene Summary: binds urotensin II and produces hypotension [RGD, Feb 2006]