

Product datasheet for **RN201274**

Gpr180 (NM_001006994) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Gpr180 (NM_001006994) Rat Untagged Clone

Tag: Tag Free

Symbol: Gpr180

Synonyms: MGC94555

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN201274 representing NM_001006994

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGACAGTGAACCACTGAACACAACCTCACAGTGTCCAGCTCCCGCGCCCCAGACGTGGCGAGTGT
TCTATGCAGACAAGTACACCTGCAAGGATGACACGGAGAGTCCCCAGGAGGAGACATCCCCTTTGAAAT
GGTGCTCCTCAACCCGGATGCCGAAGGAATCCACTGGACATTTTAGCGCCAGAGAATCTGGGCTCCAC
GAGTTCTTTTCTCCTAGTCTAGTGTACTTTGTGACTGCGTGCATCTATGCGCAGTCCCTGTGGCAGG
CTATGAGGAAAGGAGGACCCATGCACACCATCTTAAAGGTCCTCACCACCGCACTGCTGCTGCAAGCTGC
TTCAGCCTTAGCTAATTACATTCATTTCTCCAGGTAATCCAGGGATGGAATAGGGCTGCCTCTCATGGGA
AGCTTGGCAGAAGTTTTGACATTGCCTCCCAAATTCAGATGCTGTACCTACTTCTGAGCCTGTGCATGG
GCTGGACAATTGTCCGCATGAAGAAGTCTCAGAGCAGACCGCTGCAGTGGGACTCCACCCCGCGTCCAC
GGGCATGGCCGCCTTCATTGTGATCACACAGAGCGTTTTGCTGCTGTGGGAGCAGTTTGAAGACACCACT
CACCACAGCTCGCACTCACACCTCAGCCTAGCCGGCGTCTTGTGATCATCTTACGGATTTGCCTAGCGC
TGTCGCTGGGCTGTGGACTCTATCAGATCATCAGATGGAGAGGAGCACACTCAAGAGAGAGTTCTACAT
CACGTTTGCCAAAGGCTGCATCCTGTGGTTCCTGTGCCAGCCAGCACTTGCCTGCATCGCTGTCTTTTC
AATGACTACCAAAGAGATAAGGTTATCACAGTAGGTGTATCCTTTGCCAGGCTGTGGCCATGGTCATTC
TGTACAGACTTTTCTGTCCACAGTCTGTACTGGGAGGTCTCCTCACTCTCCTCAGTAACACTACCCCT
CACCATCTCGTCTGCACAAAGGGCGCCCCACTTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001006994
Insert Size:	1020 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001006994.1, NP_001006995.1</u>
RefSeq Size:	1835 bp
RefSeq ORF:	1020 bp
Locus ID:	306165
Cytogenetics:	15q24