

Product datasheet for **RN201286**

Car7 (NM_001106165) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Car7 (NM_001106165) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Car7
Synonyms:	Ca7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201286 representing NM_001106165 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGACCGGCCACCACTGTTGGGGCTACGGCCAGGACGACGGCCCTTCAAATTGGCACAAGCTGTATCCCA
TCGCCCAGGGAGACCGCCAGTCACCCATCAATATCATATCCAGCCAGGCAGTGTACTACCCAGCCTACA
GCCACTGGAGCTTTTCTATGAAGCCTGCATGTCCCTCAGCATACCAACAATGGCCACTCTGTTCAAGT
GACTTCAATGACAGTGATGACAGAACCCTGGTTGCCGGGGCCCTCTGGAAGGGCCCTATCGTCTCAAGC
AGCTCCACTTCCACTGGGCAAGAAGCGAGATGTGGGCTCAGAGCACACAGTGGATGGCAAGTCCTTCCC
CAGCGAGCTTCATCTGGTTCAGTGAATGCCAAGAAGTACAGCACTTTTGGGGAGGCGGCTGCAGCCCT
GATGGCCTGGCTGTGGTTGGCATCTTCTGGAGACAGGAGACGAGCACCCGAGCATGAACCGCCTGACAG
ATGCACTCTACATGGTTCGGTTTAAGGACACCAAGCCAGTTTCAGCTGCTTCAACCCCAAGTGCCTGCT
GCCCACCAGCCGGCACTACTGGACCTATCCTGGCTCCCTAACCACACCCCACTCAGCGAAAGTGTCACT
TGGATTGTGCTCCGAGAGCCTATCAGAACTCTGAGAGGCAGATGGAGAAATTCGGAGCCTGCTTTTCA
CCTCAGAGGACGATGAGAGGATTCACATGGTGAACAACTCCGGCCACCACAGCCGCTGAAGGGCCGAGT
GGTTAAAGCATCCTTCCAGTCC**TGA**

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001106165
Insert Size:	795 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106165.1, NP_001099635.1</u>
RefSeq Size:	1534 bp
RefSeq ORF:	795 bp
Locus ID:	291819
Cytogenetics:	19p14