

Product datasheet for **RN201345**

Fam3c (NM_198771) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam3c (NM_198771) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fam3c
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN201345 representing NM_198771 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGGTAGCAGGAGCCGAAAGTTGGTGGTAGCCGTGGCAGTGTCTTACTGACCTTCTATGTTATTT
CTCAAGTTTTCGAAATTTAAATGGATGCAAGTTTAGGAAGTCTATTTGCCGATCTGCGTTGGACTCAGC
CATTTCGCTCAACGAAGCCTCCAGGTACAAATGTGGGATCTCAAAGCGTGCCCGGAGAAGCATTGCT
TTTAAATGGCGAGCGGAGCAGCCAATGTAGTAGACCCAAGATCTGCCTGGAAGATAATGTTTTGATGA
GCGGTGTGAAGAACAATGTTGGGAGAGGGATAAATGTTGCCCTGGTGAACGGGAAAACGGGAGATGTGAT
AGATACCAAATATTTTGACATGTGGGGAGGAGATGTGGCACCATTATTGAGTTTTGAAGACCATAACAG
GACGGGACAGTAGTGTTAATGGCAACGTACGATGATGGAGCAACCAAACCTCACGGAAGAGGCACGGCGAC
TCATCGCCGAACGTTGGGAGCACATCCATCACCAGTCTTGGCTCCGAGATAACTGGGTCTTCTGTGGTGG
CAAGGGCATTAAAGACAAAGAGCCCCTTCGAACAGCACATAAAGAACAATAAGACACAAACAAGTACGAG
GGGTGGCCCGAGGTGGTCGAGATGGAAGGATGCATCCCCAGAAGCAAGACTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_198771
Insert Size:	684 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198771.1 , NP_942066.1
RefSeq Size:	1367 bp
RefSeq ORF:	684 bp
Locus ID:	312159
UniProt ID:	Q810F4
Cytogenetics:	4q22
Gene Summary:	human homolog is a member of a cytokine-like gene family [RGD, Feb 2006]