

## Product datasheet for **RN201446**

### Mapre1 (NM\_138509) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mapre1 (NM\_138509) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Mapre1  
**Synonyms:** Eb1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN201446 representing NM\_138509  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGTGAATGTGTACTCGACGTCAGTCACCAGTGATAACCTCAGTCGACATGACATGCTGGCTTGA  
TCAATGAATCTCTGCAGTTGAATCTGACAAAGATAGAACAGCTATGTTCCAGGGGCTGCATACTGTCAGTT  
CATGGACATGCTCTTCCCGGCTCCATTGCCTTGAAGAAAGTCAAATCCAGGCTAAGCTAGAACATGAA  
TATATCCAGAACTCAAATACTACAAGCCGGCTCAAGAGGATGGCGTTGACAAAATAATTCCCTGTGG  
ATAAATTAGTAAAAGGAAAATTTCAAGACAATTTGAATTCGTTCAATGGTTCAAGAAAGTTTTTGTATGC  
AAATTATGATGAAAAGAATATGATCCTGTCGCTGCCAGACAAGGTCAAGAAACGGCAGTGGCTCCTTCT  
CTTGTCGCTCCAGCTTTGAGTAAACCGAAGAACTCTCGGCTCCGGTAGTGCAGCCCCACAGAGACCCA  
TTGCAACACAGAGGACTACTGCAGCTCCTAAGGCTGGCCCCGGGATGGTGCAGAAAGAAATCCTGGCATGGG  
CAATGGAGATGACGAAGCTGCTGAATTGATGCAGCAGGTCAAAGTATTGAAGCTTACTGTTGAAGACTTG  
GAGAAGGAGAGAGATTTCTACTTCGGCAAGCTAAGGAACATTGAATTGATTTGCCAGGAGAATGAGGGGG  
AAAATGACCCTGTACTGCAGAGGATTGTAGATTTCTTTATGCCACAGATGAAGGCTTTGTGATACCTGA  
CGAGGGGGGCCACAGGAGGAACAAGAAGAGATTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_138509  
**Insert Size:** 807 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_138509.3</a> , <a href="#">NP_612518.2</a>
<b>RefSeq Size:</b>	2310 bp
<b>RefSeq ORF:</b>	807 bp
<b>Locus ID:</b>	114764
<b>UniProt ID:</b>	<a href="#">Q66HR2</a>
<b>Cytogenetics:</b>	3q41
<b>Gene Summary:</b>	expressed specifically at growing microtubule distal tips [RGD, Feb 2006]