

Product datasheet for RN201479

Spryd4 (NM_001037765) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spryd4 (NM_001037765) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Spryd4
Synonyms:	RGD1306649
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201479 representing NM_001037765 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCGCTGCCCTTTGCACGTTGTTGGAAGTTGTACCGCTGGGAACCAAGCGATGGGGGTACCGCCG
 GGGAAATCCCGAGAGGCATCAGTTTCAAAGTGAAGAAAAACGGCCACAGCAGCCTGGCACTCTTCAG
 AGGTGACACGGGTGTCAAATATGGCTTGGTGGGACTGGAACCCACCAAGGTGGCCCTGAAGTTGGAGCGC
 TTCCGGGACTGGGCAGTGGTGTCTCGACACACCACAGTCACCAAGTGGCAGACACTACTGGGAGGTAAGT
 TGAAGCGCTCCCAGCAGTTCCGGATAGGAGTGGCAGATGTCGACATGTCCGGGACAGCTGCGTTGGTGC
 TGATGATCGTTTCTGGGTGTTCTCTATGCCAGCGCAAGTGGCAGCAGATGTTGGCAAACGAAAAAGCC
 CCAATTAAGGGTATCGGGCAACCAGAGAAAGTGGGGCTGCTGCTGGACTATGAGGCAAAGAAGCTGAGCC
 TAGTGGATGTGAATCGGATCTCTGTCTATCCACACACTACAGACAGATTTCCGGGGCCAGTGGCACCTGC
 GTTTGCTCTTTGGGATGGAGAACTGCTAACGCACTCAGGACTTGAAGTGCCCAAGGGCCTC**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001037765
Insert Size:	624 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037765.1, NP_001032854.1</u>
RefSeq Size:	1518 bp
RefSeq ORF:	624 bp
Locus ID:	288772
UniProt ID:	<u>Q4FZT8</u>
Cytogenetics:	7q11