

Product datasheet for **RN201490**

Mrnip (NM_001109030) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mrnip (NM_001109030) Rat Untagged Clone
Tag: Tag Free
Symbol: Mrnip
Synonyms: RGD1564036
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN201490 representing NM_001109030
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGCCGGCTCAGCAGTCTCGGGTTCTACGGTGTGCAGCTGCCACATCTTTCAGGCACATCAGGTGA
AGAAGAGTCTCAAGTGGACATGCAAGGCTTGTGGAGAGAAGCAGTCTTTCGTGCGGGCTTATGGCGAAGG
CTCTGGTGTGACTGTAGACGCCATGTCCAAAAGCTAAATCTGTCTCAAGGACAAGCGTCAAGAGCTGTCA
CTCAGGTCTGTGGAGGAAGCTGTAATGCCAGTGAAGAGGAAAATGCAGGCCCTCAGCAGCTAAGGATG
GAAGTCAGCAGGCATTGTCTGAGCCCCTCGAGTCCCGTTGGCTGAAGTATCTGGACCAGGACTGTGAGGG
CCAGGAGCTGGACAGAGGAAGACCAGCTCTTAAGACACAGCTCTCAGCCAGTGCTAAGTGGCCAAATCTG
CCCTGCAGTCCAGCCGACCCAGGAAGAGGAAATGGGATCAGAGCACAGGGCAGCCTGCACACAGCCCGC
GTGTCCAGGGTGTGGACAGTGGTGAAGCAACTTTCAACACCAGGACTCCACTGGCCTGTCTGGGACAGA
GCAGCAAGGCATTAGCCCTGCCCTCAGTACCACAAACCACATCAGGGAGCTCGGTTTTCTAGGTGGAAA
CTACCAAGTCTGTACGCAAGTTAATGCCCATCTTCTAAATGGGCACGGTTTTCTCCTATCGCCTGAGA
ACAGTCCACAAGTGGATGGGAAACCAACAAAGTCTTTGCAGGAGGGCACCAGGCGATCTGGCCTAGCACA
AATTGAGCAAGGGACTCTTGAGACCAAGACTCCAGGAGAAGGCCACTTCATCAGAGACCCTGCCACCATC
CGGCCGCCTTAACCACTCACACTACCCTCCACAGCCAGACAAGCCTGACAGGAAGACCCAGAGCAGC
CATGGGCCACGGGACTCCTCAAGCAGATGGTGGGCTTTTGGCTCAAGGGACACAGAAGGCTCCATCCAT
GCAACTGCACAACCTCTTTACAACCGGAGAGGACTTTGAGGATGACCTGTAG

ACGGTACGGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001109030



[View online »](#)

Insert Size:	1032 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109030.2</u> , <u>NP_001102500.2</u>
RefSeq Size:	1173 bp
RefSeq ORF:	1032 bp
Locus ID:	497895
Cytogenetics:	10q22