

Product datasheet for **RN201557**

Trim52 (NM_001106056) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Trim52 (NM_001106056) Rat Untagged Clone
Tag: Tag Free
Symbol: Trim52
Synonyms: RGD1565615
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN201557 representing NM_001106056

Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCCGCCCTACAGGGCCTCCAGCCCTATGCAGTCACTTCGGGAAGAAGCAGTATGTGCCATCTGTC
 TGGATTACTTTAAGGACCCCGTCCATTGGCTGTGGGCACAACCTTGCCGAGGGTGTGTGACCCAGCT
 GTGGGGCAAGGAAGATGAACAGGACCGGGAACCGCCTGCTGTAGGAACACGTCATTAGAGAGGTTTTG
 TTTCATAGGTACACAGAACAGGAGGTGTTTCAGCACCGAGCACATGCTGGACGCTGGGCTGCTCATAGCC
 ATAGAAGGCATCGGGCAATGCAGACTCTGTGTGGGATGATGAGGAAGAAGACTGGAACAGTTTACAAGG
 ACTGGTACATGACCTGAGAATTAGGGTTTTTCCAGAAGAAAGAGATGAACACCCCCACGATGGCCACCAA
 TACCATCACTATGGCCGATACCGACACCGTCCAGTTTTCCGCCGTGGTCCCCACATCCACCTGTGCGTC
 GGCAGTCTATCCAGATGCCCGACACGATCTCCACCTCGAGTACGATCTCCACCTCGAGTACATTCTCC
 ACCTCGAGTACGTTCTCCACCTCGTCCACATCACAGGTCTTCAGCTGCCCAATGCCGAAGGACATTT
 CCAAGTCGCAGTTTTCGACCAATTTGCAGCTGGCCAACATGGTCCATATAATTGCCAGATTGGCCATA
 CGCCATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001106056
Insert Size: 708 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106056.1, NP_001099526.1</u>
RefSeq Size:	880 bp
RefSeq ORF:	708 bp
Locus ID:	290458
Cytogenetics:	15q22