

Product datasheet for **RN201574**

Mtap (NM_001047867) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mtap (NM_001047867) Rat Untagged Clone
Tag: Tag Free
Symbol: Mtap
Synonyms: Cc1-6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN201574 representing NM_001047867
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAAGTCGGCGGCAGATGCTGGCTGATACTCTGGCTTTGGGATTGATTGCCCGCAGTCTCCCCTTG
GTCCAGCGCTGCGGCCCGCGTGGACTTCCCCGGTGCAGACATGGCCTCGGGCTCCGCTGCACGGCAGT
GAAGATTGGAATAATTGGTGAACAGGCTTGGATGACCCTGAAATTTAGAAAGGAAGAACAAGAAAAATAT
GTGGATACTCCATTTGGCAAGCCATCTGATGCCTAATTTGGGGAAGATAAAGAACGTCGACTGTGTGC
TTCTTGCCAGCATGTACGTTTTAGATCGTATGTGTGAGTTCCTGCTGCTGCCTTTCTGTGTCCAGCTTAA
TTCACCTCACGAGATTCCTCCAGGCTCATCTTGTACACCTCTTATCTGTTTGCCTTTGTTGCCCGCTGG
GAAAGGCGTGTCCGCTGGTTAGCGAGAGTTCCTGTGAGCGGGCCTTTGTGCCGAGCTCAGCATCTACTT
CCGTTCCAGTGGAAAGAAACCAACATTCGAAGGAGATAACGCAAAGATCCAAGCATCCAGGTTCTCCAC
GACAGTGCAACAGCAGGCTATCCGAGAGGAGTGCCAGACATCCCTGAGGCCTCAGACCTTCTATGATGGA
AGTCATTGCAGTGCCAGAGGAGTGTCCACATCCCAATGGTGAACCATTTGCCCAAAACAAGAGAGG
TCCTCATAGAAATGGCTAAGAAGCTTGGACTCCGGTGCCATTCAAAGGGGACGATAGTCACAATTGAGGG
GCCTCGGTTTAGCTCCCGGCAGAAAGCTTTATATTCCGTAAGTGGGCGCAGATGTTATCAATATGACC
ACTGTTCCAGAGGTGGTCTTGCAGAGGAGGCGGATTGTGCTACGCAAGCATTGCCATGCAACCGATT
ATGATTGCTGGAAGGAGCATGAAGAAGCAGTGTGAGTGGATGGGTTTTAAAGACCATGAAAGAAAATGC
CAATAAGCCAAAAGTTTACTCCTAACTACTATACTCAGATAGGGTCAATGGAATGGTCAGAAACGCTC
CGTAACCTAAAGGTGGTTGCTGTCAAGCTGTATGGAGAAAAAAGAGAGAGCTGCACCAAAAATCGTCTG
GCGGTATGCTGCAGATGCTTGCCTTCCAAGGACTCTGGACAATTGGCAGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_001047867
Insert Size:	1176 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001047867.1, NP_001041332.1</u>
RefSeq Size:	1182 bp
RefSeq ORF:	1176 bp
Locus ID:	298227
Cytogenetics:	5q32