

Product datasheet for **RN201580**

Bdkrb1 (NM_030851) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Bdkrb1 (NM_030851) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Bdkrb1
Synonyms:	b1bkr; Bdkrb; BKR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201580 representing NM_030851 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGTCCGAGGTCTTGTGGAGCTCCAGCCCTCTAACCGAAGCCTGCAGGCCCTGCCAACATTACCT
CCTGCGAGAGTCCCTAGAAGACTGGGACCTGCTGTATCGGGTGTGCCAGGGTTCGTCACACTATCTG
CTTCTTTGGCCTCTGGGAACCTTTAGTCTTGTCTTCTTCTTTGCCTTGGCGACAGTGGTGGTGG
CAGCAACGACAGAGGCAGCAGCGCTTAACCATAGCGGAAATCTACCTGGTAACCTGGCGGCTCCGATC
TGGTGTGTCTGGCCTGCCCTCTGGGCAGAGAACATCGGGAACCGTTTCAACTGGCCCTTCGGAAC
TGACCTCTGCCGGTGGTCAGCGGGTCAATCAAGCCAACTGTTTGTACGATCTTCTGGTGGTGGCT
ATCAGTCAGGACCGCTACAGGCTCCTGGTATACCCATGACCAGCTGGGGTACCGCGCGACGGCAAG
CCCAAGCTACGTGCCTGCATCTGGGTAGCCGGGGTCTCTTGAGCATCCCCACATTCTTCTACGCTC
TGTTAAAGTCGTCCCGATCTGAACGTCTGCCTGCATCCTGCTTTCCCCACGAGGCTTGGCACTTT
GCAAGGATGGTGGAGTTGAACGTTTTGGGTTCTCTCTCCCGTACTGCTATCATCTTCTCAACTATC
ACATCTGGCCTCCCTGAGAGGACAGAAGGAGGCCAGCAGGACTAGGTGTGGGGTCCCAAGGGCAGCAA
GACGACGGGGCTGATCCTCACACTGGTAGCCTCCTTCTGGTCTGCTGGTGCCTTACCCTTCTTCGCT
TTCCTGGATTTCTGGTCCAGGTAAGAGTGATCCAGGACTGCTCCTGGAAGGAGATCACAGACCTGGGCC
TGCAGCTCGCAAATCTTTGCCTTTGTCAACAGCTGCTTGAACCACTGATTTATGTCTTCGACAGGCC
GCTCCTTAAGACCAGGGTCTGGGACTTTA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1981_g02.zip



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_030851
Insert Size:	1014 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_030851.1 , NP_110478.1
RefSeq Size:	1312 bp
RefSeq ORF:	1014 bp
Locus ID:	81509
UniProt ID:	P97583
Cytogenetics:	6q32
Gene Summary:	may play a role in regulation of inflammatory response [RGD, Feb 2006]