

## Product datasheet for **RN201795**

### **Ccr7 (NM\_199489) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ccr7 (NM\_199489) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ccr7  
**Synonyms:** MGC108519  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN201795 representing NM\_199489  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGACCTGGGAAGCCACGAAAAACGTGCTGGTGGTGGCTCTCCTGGTCATTTCCAGGTGTGCTTCT  
GCCAAGATGAGGTCACAGACGACTACATCGGCGAGAACCACCCTGGACTACACCCTGTATGAGTCGGT  
GTGCTTCAAGAAGGATGTGCGGAACCTTAAGGCCTGGTTCCTCCCTCTCATGTACTCAGTCAATTTGCTTC  
GTGGCCCTGCTAGGCAATGGGCTGGTGGTGTGACATACATCTATTTCAAGAGACTGAAGACCATGACGG  
ATACCTACCTGCTCAACCTGGCCGTGGCAGACATCCTCTTCTCATGATCCTTCCCTTCTGGGCCTACAG  
CGAAGCCAAGTCCTGGATCTTTGGTGCCTACCTGTGTAAGAGCATCTTTGGCATCTACAAGTTAAGCTTC  
TTCAGCGGGATGTTGCTGCTCCTGTGTATCAGCATTGACCGCTATGTGGCCATCGTCCAGGCCGTGTCAG  
CCCACCGGCACCGCCCGCGTCTTCTCATCAGCAAGCTGTCTGTATAGGCATCTGGACGCTGGCCTT  
TTTCTTTCTATCCCTGAGCTGCTCTACAGCGGCTCCAGAAGAAGCAGCGCGGAGGACACGTGGAGATGC  
TCCCTGGTCAGTGCCCAAGTGGAGGCCTTGATCGCCATCCAAGTGGCCAGATGGTTGTTGGGTTGTAC  
TGCTATGCTGGCTATGAGTTTCTGCTACCTGGTTATCATCCGCACTCTGCTCCAGCGCGAAACTCGA  
GCGGAACAAGGCCATCAAGGTGATCATCGCTGTGGTTCGTAGTGTTCGTCGTCTTCCAGCTGCCCTACAAT  
GGGGTGGTCTGCCCCAGACCGTGGCCAATTTCAACATCACCAATAGCAGCTGCGAAGCCAGCAAGCAGC  
TCAACATTGCCTATGACGTACCTACAGCCTGGCCTCCGTCGCTGTGTCAACCCTTCTTGTATGC  
CTTCATCGGCGTCAAGTTCGCGACGACCTCTTCAAGCTCTTCAAGGACTTGGGCTGCCTCAGCCAGGAA  
CGGCTCCGGCAGTGGTCTTCTGCCCATGTACGGCACACGTCCTGAGCATGGAGGCGGAGACTACCA  
CCACCTTCTCCCC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_199489
<b>Insert Size:</b>	1137 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_199489.4, NP_955783.1</u>
<b>RefSeq Size:</b>	1215 bp
<b>RefSeq ORF:</b>	1137 bp
<b>Locus ID:</b>	287673
<b>Cytogenetics:</b>	10q31
<b>Gene Summary:</b>	mouse homolog is a chemokine receptor; involved in the immune response to viruses [RGD, Feb 2006]