

Product datasheet for **RN201800**

Fcgr3a (NM_207603) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fcgr3a (NM_207603) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fcgr3a
Synonyms:	Fcgr4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201800 representing NM_207603 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGTACCTACTACTACCAACGGCACTGCTACTTACGGTTTCCTCTGGCGTTGGAGCTGGTCTCCAAA
AGGCTGTGGTGAACCTAGACCCTGAATGGGTCAAGGTTCTCGAGGAAGACTGTGTGATCCTCAGGTGCCA
GGGCACCTTCTCCCCGAGGACAATTCTACAAATGGTCCATAACAAAAGCCTCATCTCGCACCAGGAC
GCCAATATGTCATCAAAGTCCAGAGTTAAGGACAGTGAATGTACAGATGCCAGACAGCCTTCTCCG
CGCTCAGTGACCCGGTGCAGCTAGACGTCCATGCAGACTGGCTATTGCTTCAGACCACTAAGCGGCTGTT
CCAGGAGGGGACCCATTTCGTCTGAGATGCCATAGCTGGCGAAACACGCCTGTATTTAAGTTACCTAT
TTACAGAAATGGCAAAGGCAAGAAGTATTTCCATAGGAATTCGAACTCAGTATTTCAAAAGCTACGCACG
CCGACAGTGGTTCCTACTTCTGCAGAGGGATCATTGGACGCAACAACATATCTTCAGCATCCTTGAGAT
AAGCATAGGAGATCCGACGTCTCCCTCCAGCTTCTACCGTGGCATCAAATCACTTTCTGCCTGCTGATA
GGACTCTGTTTGCAATAGACACAGTGTGATTTCTCGGTGCAGAGGAGTCTTCAAAGTTCCGTGGCAG
TCTATGAGGAACCCAACTTCACTGGAGCAAGGAACCTCAGGACAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms:	https://cdn.origene.com/chromatograms/ja1334_g10.zip
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_207603
Insert Size:	750 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207603.1</u> , <u>NP_997486.1</u>
RefSeq Size:	1178 bp
RefSeq ORF:	750 bp
Locus ID:	304966
Cytogenetics:	13q24