

Product datasheet for **RN201905**

Spsb4 (NM_001106849) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spsb4 (NM_001106849) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Spsb4
Synonyms:	RGD1562307
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201905 representing NM_001106849 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGGTCAGAAGCTCTCTGGGAGCCTCAAGTCTGTGGAGGTGCGAGAGCCAGCGCTTCGGCCAGCCAAGC
 GGGAGCTCCGAGGCTTAGAGCCTGGGCGGCCAGCACGGCTAGACCAGCTGCTAGACATGCCTGCAGCGGG
 CTTGGCTGTACAGCTGAGACACGCGTGAACCTGAGGACCGATCCCTCAATGTCTTCGTAAGGATGAC
 GACCGCCTCACCTCCACCGACACCCGGTAGCCAGAGCACAGATGGCATCCGAGGCAAGGTTGGCCATG
 CCCGGGTCTGCATGCATGGCAGATCCATTGGCCAGCTAGGCAGCGGGGTACCCACGCTGTGGTGGGCGT
 GGCCACAGCCCGTGCCCTTTGCACTCTGTGGTTACACCGCACTGGTAGGCAGTGACTCCGAGTCGTGG
 GGCTGGGACCTGGGTCGTAGCCGCCTCTACCACGACGGCAAGAACCGGCCTGGGGTGGCATAACCTGCCT
 TCCTGGGGCCAGATGAGGCCTTTGCTCTGCCGATTCTTATTGGTTGTATTGGACATGGATGAGGGCAC
 GCTCAGCTTCATCGTGGATGGCCAGTACTTGGGCGTGGCATTCCGTGGCCTCAAGGGCAAGAAATATAC
 CCGGTGGTGAGTGCTGTATGGGGTCACTGCGAAGTCACCATGCGCTACATCAATGGCCTTGACCTGAGC
 CTCTGCCGTTGATGGACCTATGCCGGCGATCCATCCGCTCAGCCCTGGGCAGACAGCGCCTGCGGGACAT
 CGGCTCCCTGCCCTGCCTCAGTCTCTTAAAACTATCTGCAGTACCA**TGA**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-RsrII
ACCN:	NM_001106849
Insert Size:	822 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106849.2, NP_001100319.1</u>
RefSeq Size:	2802 bp
RefSeq ORF:	822 bp
Locus ID:	300950
Cytogenetics:	8q31