

Product datasheet for RN201930

RT1-Ha (NM_001008848) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RT1-Ha (NM_001008848) Rat Untagged Clone
Tag:	Tag Free
Symbol:	RT1-Ha
Synonyms:	RT1.Ha
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201930 representing NM_001008848 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTCCAGACTAAAGCTGAGGCCCTGAGAGCCTCCTTCTGGCCTTCTGCTGAGTCTCCAGGATCTG
GGGCCATCAAAGCGGAGCACGTAGCAACGTACGACATATTTGTGCAGACGCAGCACCCGCTGGTGAATG
CATGTACGCGTTTGACGGTGATGAGCAGTTCTACGTGGATCTTGACAGGGAGGAGACGGTCTGGGACCTG
CCGGAGTTCAGCAAGGTGATGCGTTTGACGCCAGGAGGCTCTCCCTACATCGACAGCCTGAAGAACA
ACTTGAGAGCCTTGATCCAGAGATACAACGTACCCAGGACCCAGCGAGCCTCCCGAGGTGACCGTGT
CCCTAAGAACCCGGTAGAGCTGGGCCAGCCAAATGTCCTCATCTGCCACATCGACAGGCTCTTCCCCCA
GTGCTCAATGTACATGGCTACGTAACGGGCAGCCGGTACCGAAGGAACCTCAGAGACCGTCTTCTGG
CCAGCACAGAACTCAGATTTACAAGTTCACCTACCTGACCTTATTCCACAGCGCAGGATGTGTATGA
CTGCAGGGTGGAGCACTGGAGCCTGGACCAGCCGAGTCTCAGTCACTGGGAAATGCAGGAGCCAGTCCAG
GTGTCGAGACAATGGAGACCAGGTCTGCGCCCTGGGCTGGTGGTGGGCTGGTGGGCTCCGTGGTTG
GTAGTGCCTCATCGTAAGGGCTCCGCATGGCGGCCATAACCCCTCGGACCCAGGAACCCCTG**TGA**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-RsrII
ACCN:	NM_001008848
Insert Size:	765 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008848.2</u> , <u>NP_001008848.1</u>
RefSeq Size:	839 bp
RefSeq ORF:	765 bp
Locus ID:	24986
Cytogenetics:	20p12
Gene Summary:	may play a role in antigen presentation [RGD, Feb 2006]