

Product datasheet for RN201978

Olr200 (NM_001000189) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr200 (NM_001000189) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Olr200
Synonyms:	ratchr1-163138220-163139185_ORF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN201978 representing NM_001000189 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGGTACTGCCTTGCATGAAACCAACAGCTCAGACTTCCACGTGACCGAGTTCATTCTGCTGGGATTCC
 CAGGCATTCATGACTCCAGCACTGGATCTCCTTGCTATGGCTCTGCTTACATTGTAGCTCTTGGTGC
 CAATATCCTTATCTTGATCACCCTCTACTGTGAGCCAAGCCTGCATCAGCCAATGTACCAGTTTCTTGGT
 ATCTTGGCTGCAGTAGATATTGGCTTGGCTACCACCAGCACCCCAAGATCCTGGCCATCTTGTGGTTTG
 ATGCCAAGACCATCAGCCTCCCTGAGTGTCTTGTCAAATCTATGCAATCCATACTTTATGTGCATAGA
 GTCAAGCATCTTCCTGTGCATGGCAATAGATAGATACGTAGCCATTTGTTATCCCTCCGGTATCCTTCC
 ATAGTCACTGAAGCTTTTGAATCAAAGCCACACTGACCATGCTGCTAAGGAACAGCCTGTTACCATGC
 CAGTCCCTGTACTAGCTGCCAGAGACAGTACTGCTCCAGAAATGAAATAGATCACTGCCTTTGCTCTAA
 CTTGGGCGTTATTAGTCTGGCCTGTGATGATATCACTGTCAACAGGTTCTACCAGTTGACCTTGGCCTGG
 ATTATGGTTGGAAGTGACATGATTCTGGTCTGTGCATCCTATGCTTTGATTATTTGCTCAGTACTGAGGC
 TGAACCTCACTGAGGCAATATCTAAGGCCCTCAGTACCTGCAGCTCCACCTCATTCTCATTGTTTGA
 CTACCCAGCTATTGTTATATTATCCGTCACTCACCTAGCAGGAAGAAGGGTTCCCTCATCCAGTTCTC
 CTCAATGTGATGCATATTGTCATCCCCCATCTCTTAACCCAGTGGTATATGCCCTTAGGACACAGGAGC
 TGAAGGTAGCCTTCCAGAAGTACTTGGCTTTAGTGAGTATATGTCCAGGAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001000189


[View online »](#)

Insert Size:	966 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000189.1, NP_001000189.1</u>
RefSeq Size:	966 bp
RefSeq ORF:	966 bp
Locus ID:	293340
Cytogenetics:	1q32
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]