

## Product datasheet for **RN202037**

### Cplx2 (NM\_053878) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cplx2 (NM_053878) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Cplx2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202037 representing NM_053878 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGACTTCGTCATGAAGCAAGCCCTCGGAGGGGCCACCAAGGACATGGGAAGATGCTGGGGGAGAGG  
AGGAGAAGGACCCAGATGCACAGAAGAAGGAGGAGGAGCGGCAGGAGGCCCTGAGGCAGCAGGAGGAAGA  
GCGCAAGGCCAAACATGCCCGCATGGAAGCGGAACGTGAGAAGGTCGGCAGCAGATTCGAGATAAGTAT  
GGGCTGAAGAAGAAGAAGAGAAGGAGGAGGAGGAGAAGGCAGCCCTGGAACAGCCCTGCGAGGGAAGCC  
TGACCAGACCCAAGAAGGCCATCCCTGCAGGCTGTGGGACGAGGAGGAGGAAGAAGAGGAGGCATCCT  
GGACACGGTGCTCAAATATCTGCCGGGGCCACTGCAGGACATGTTCAAGAAG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_053878
Insert Size:	405 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_053878.1](#), [NP\\_446330.1](#)

**RefSeq Size:** 912 bp

**RefSeq ORF:** 405 bp

**Locus ID:** 116657

**UniProt ID:** [P84087](#)

**Cytogenetics:** 17p14

**Gene Summary:** binds the assembled SNARE (soluble N-ethylmaleimide-sensitive factor attached protein receptor) complex; modulates synaptic vesicle exocytosis [RGD, Feb 2006]