

Product datasheet for **RN202065**

Tmem158 (NM_057212) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Tmem158 (NM_057212) Rat Untagged Clone

Tag: Tag Free

Symbol: Tmem158

Synonyms: Ris1

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN202065 representing NM_057212
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGCTACCCCTTCTCGCCGCGCTACTGGCCGCCGCTGCCGCTGCCGCCGGCGCGCGGGGCCACGG
ATGCGCCCGGCTTTCTGGGACGCCCCAAATGCTTCGGCCAATGCTTCCTTCACCGCGAACACTCCAC
CCCGCTCTGCTGGCATCGGCGGCATCCGACCCCCAGAGCGCGCGGGCCAGAGGAGGCTCCTGCGGCG
CCCTGCAACATCAGCGTGCAACGGCAGATGCTGAGCTCCCTGCTCGTGCCTGGGCGCGCCGCGGGCT
TGCACTGCGACCTGCTGCTCTTCTCTACCAACGCGCACGCGCGCCTTCTTTGCCGCGCCTTCCACCG
CGTCGGGCGCGCGCTGCTCATCGAGCATCTGGGGCTGGCGGCCGGCGCGCGCAGCAAGACCTGCGCCTG
TGCGTGGGCTGTGGCTGGGTGCGCGGCCGCTGCGAGCTCCCGCGGGGCTCCACCGCTCTGCCTGCCT
ACCCCGCGGCCGAACCGGGGCGCTGTGGCTGCAGGGCGAGCCTCGTCACTTCTGCTGCCTGGACTTCAG
CCTGGAGGAGCTGCAGGGCGAGCCGGGCTGGCGGCTGAACCGCAAGCCCATCGAGTCCACGCTGGTGGCC
TGCTTCATGACCTTGTCTCATCGTGGTGTGGAGCGTGGCCGCCCTCATCTGGCCGGTGCCCATATTGCCG
GCTTCTTGCCCAACGGCATGGAGCAGCGCCGACCAACCGCGCGCCCCGCGGCCGCTCCTGCAGCCGT
ACCCGCGGGGACTACCGCCGAGCCGCCGCGCTGCTGCAGCAGCTGCCGCGCCGCGGCCGTACCTCG
GGGTGGCTCCAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_057212

Insert Size: 858 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_057212.1, NP_476560.1</u>
RefSeq Size:	999 bp
RefSeq ORF:	858 bp
Locus ID:	117582
UniProt ID:	<u>Q91XV7</u>
Cytogenetics:	8q32
Gene Summary:	receptor that binds brain injury-derived neurotrophic peptide; abolishes the protective effect of brain injury-derived neurotrophic peptide [RGD, Feb 2006]