

Product datasheet for RN202112

Pctp (NM_017225) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pctp (NM_017225) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pctp
Synonyms:	PC-TP; stARD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202112 representing NM_017225 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCGGGGCCCGCAGCCCACTTCTCGGACGAGCAGTTCGGGAGGCTTGTGCCGAGCTCCAGAAGCCGG
CGCTAACTGGGGCGGATTGGCAGCTCCTGGTAGAAGCCTCAGGCATAACCATCTACCGGCTGCTGGACCA
GTCAACTGGACTTTACGAGTATAAAGTGTGGAGTCTGGAAAGTTGCATACCATCTCTACTCGCAGAC
GTTTACATGGACTTAGACTACAGGAAAAATGGGACCAGTACGTGAAAGAACTCTATGAGAAGTCGTTTG
ACGGACAGATGGTGGCGTACTGGGAAGTGAAGTACCCTTCCCTCTGTCCAACAGAGATTATGTCTACAC
CCGCCAGCGGCGAGACCTGGATGTGGACGGGAGGAAGATCTACGTGGTCTGGCCAGAATATCTCTGTA
CCTCAGTTTCCGGAGAAGTCTGGGGTAATCCGAGTGAAGCAGTACAAGCAGAGCCTGGCGATCGAGAGCG
ATGGCAAGAAGGGGAGCAGAGTTTTATGTACTACTTTGATAACCCGGGTGGCCAAATCCGTCCTGGCT
CATTAAGTGGGCAGCAAGAATGGAGTTCCAAGCTTTTGAAGACATGGTGAAAGCGTGTCTCAGAACTAC
CACAAGAAAAC**CTAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_017225
Insert Size:	645 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017225.3</u> , <u>NP_058921.1</u>
RefSeq Size:	1771 bp
RefSeq ORF:	645 bp
Locus ID:	29510
UniProt ID:	<u>P53809</u>
Cytogenetics:	10q26
Gene Summary:	plays a role in intermembrane transfer of phosphatidylcholines [RGD, Feb 2006]