

Product datasheet for RN202141

Trib1 (NM_023985) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trib1 (NM_023985) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Trib1
Synonyms:	Gig2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202141 representing NM_023985 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGGGTCGGTCCCGTGCCTTTGCCCTGAGCGGCGCCTCGCAGCCCCGCGGACCGGCTCTACTGTTCC
CGGCTGCCCGGGCACCCTCGCCAAACGCCTGCTGGACACGGACGACGCGGCGGCCGTGGCGCCAAAGTG
TCCGCGCCTCTGAGTGTGGAATCCCCTGACTACCTCAGTCCCCCGGCTCGCCCTGCAGTCTCAG
CCTCCACCCTCGGCGCAGGGACCGGCGGAGCTGTGTGAGCGCACCCGGGCCAGCCGAATCGCCGACT
ACCTGCTGCTACCCTAGCTGAGCGGAGCATGTGTCCCGGGCACTGTGCATCCACACCGGCCGCGAGCT
GCGCTGCAAGGTGTTCCATTAAACACTACCAAGACAAAATCAGGCCTTACACCCAGCTGCCGTCCCAT
AGAAACATTACTGGCATCGTGGAAGTGATCCTCGGGGAGAGCAAGGCCTACGTTTTCTTTGAGAAGGACT
TTGGAGACATGCACTCCTATGTACGAAGCCGAAAAGGCTTCGGGAAGAGGAGGCCGCTCGGCTCTTCAA
GCAAATGTTTCCGCGCTTGCCCACTGCCACCACTGAGCCATTGTGCTGGGGACCTGAAGCTTAGGAAA
TTCGTCTTCTCCACAGAGGAGAGGCCAGCTTAGACTGGAAAGCCTAGAAGATACGCACATGATAAAGG
GTGAAGATGACGATTGTGAGATAAACATGGCTGCCCGGCTATGTGAGCCCTGAGATCCTCAACTAC
TGGGACCTACTCCGAAAGGCAGCGGATGTTGGAGCCTTGGGGTGATGCTCTACACCCTTTTGGTCGGA
CGATACCCCTTTCATGACTCAGACCCGAGCGCCCTTTTCTCCAAAATCCGTCGGGGACAGTTTTGCATT
CCGAGCATGTTTCCCAAGCCAGGTGCCTCATCCGCAAGCCTCCTGAGACGAGAACCTTCTGAGAGACT
CACAGCTCCTGAGATCTTACTCCATCCCTGGTTTGAATACGTCTTGAACCGGGGTATGTTGACTCAGAA
GTAGGAATTCTGACCAGATTGTTCCCGAGTACCAGGACGACAATGACGTTAGTCTTTTTCTGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_023985
Insert Size:	1119 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_023985.1 , NP_076475.1
RefSeq Size:	1972 bp
RefSeq ORF:	1119 bp
Locus ID:	78969
Cytogenetics:	7q33
Gene Summary:	a nuclear factor; gene expression induced by m1-acetylcholine receptor [RGD, Feb 2006]