

Product datasheet for **RN202233**

Mrps9 (NM_001100549) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrps9 (NM_001100549) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Mrps9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202233 representing NM_001100549 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGCCCCCTGTATGTCCTGCGGCCGTGGGCTTTCCTTCGGCTTGCTCACGCGTTAGGGCTAACC
TCTGCCAGAGGCCCGGTACTGGACAGCCCCGGCGGTTGGCTGGCAGACGGGTCCAGGTCCCAGCTTTT
AAAGCTAATACATACCACAGTTGTAACAACAAAGAAAAATGCCAAGCCTCAAGACAGGAATCGTATACA
GAGGATTTTATTAGAAAGCAGATTGAAGAATTCACCTAGGGAAGAGGCATTGGCCAACATGATGGGAG
AAGACCCAGAAACCTTTACCGAAGAGGACATTGACAGAGCTATTGCCTACCTTTTCCCAAGTGTTTATT
TGAGAAGCGAGCCAGGCCAATGATGAAGCATCCTGAGCAAATTTTCCAAAGCAAAGAGCAATCCAGTGG
GGAGAAGACGGCCGTCCATTTCTTTCTTCTATACTGGCAAACAGTCGTACTACTCGCTAATGCACG
ACGTGTATGAAAAGTTCTCCAGCTAGAAAAGCACAGAGGTCCGCTTTCAGAAAACACTGAGTCCAGAGA
CCTGATTGGCAGCAGATGGCTGATTAAGGAGGAAGTGAAGAAATGTTAGTGAAAAACTGTCAGATCAA
GATTACACGCAGTTCATCCGACTGCTGAAAAAGTTACTGACATTGCCGTGCGGTCTGCCGAAGAGGACT
TTGTGCAGAGGTTCCGCCGAGTGTAACGTTCATCCAAAAGCAGCTGATCGAACCCGTTTCAGTATGA
TGAGCAGGGCATGGCCTTCAGCACCAGCAAGGCCGAGGAAGAGCGCCACAGCACGGGTGTTGTTTAT
CAGCACGGAAGCGGAAGATCCACGTGAACGGAGTAGATTATCTGCTTTACTTTCCCATCTTACAGGACA
GAGAACAGTTGATGTTCCCTTTCCACTTCTGGACCGGCTAGAAAAGCACGATGTTACCTGCACAGTCTC
CGGGGGAGGGAGATCAGCGCAGGCAGGACCATACGCTTGGCAATGGCCAAAGCTTTATGCAGCTTCATC
ACTGAGGACGAGGTTGAGTGGATGAGACAAGCCGACTACTTACTCCTGATCCACGGATCAGAGAACGGA
AGAAGCCAGGCCAGGAGGTTGCCGAAGGAAGTTACCTGGAAGAAGCGC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001100549
Insert Size:	1173 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100549.1, NP_001094019.1</u>
RefSeq Size:	1404 bp
RefSeq ORF:	1173 bp
Locus ID:	301371
Cytogenetics:	9q22