

## Product datasheet for **RN202236**

### **Cstf2 (NM\_001131014) Rat Untagged Clone**

#### **Product data:**

|                    |   |
|--------------------|---|
| Product Type:      | Expression Plasmids                     |
| Product Name:      | Cstf2 (NM_001131014) Rat Untagged Clone |
| Tag:               | Tag Free                                |
| Symbol:            | Cstf2                                   |
| Synonyms:          | MGC188397                               |
| Vector:            | pCMV6-Entry (PS100001)                  |
| E. coli Selection: | Kanamycin (25 ug/mL)                    |
| Cell Selection:    | Neomycin                                |



[View online »](#)

**Fully Sequenced ORF:** >RN202236 representing NM\_001131014  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCGGGTTGCTGTAGAGACCTGCGGTGGATCGTTCTCTGCGTTCTGTATTTCGTGGGAACATTC  
CTTATGAAGCTACTGAAGAACAACCTGAAGGACATCTTCTCTGAGGTTGGCCTGTTGTTAGTTTCAGATT  
GGTATACGATAGGGAGACAGAAAACCAAGGGTTATGGCTTCTGCGAATACCAAGACCAAGAGACAGCA  
CTTAGTGCCATGCGGAACCTGAATGGGCGTGAATTCAGTGGGAGAGCACTTCGAGTGGACAATGCTGCCA  
GTGAAAAGAACAAGAGGAACTGAAGAGCCTTGGCACTGGTGCACCTGTCATAGAGTACCTTATGGAGA  
AAGCATCAGCCCTGAGGATGCCCTGAATCCATTAGCAAAGCTGTTGCCAGTCTTCCGCCAGAACAGATG  
TTTGAAGTATGAAACAATGAAGCTCTGTGTCCAGAATAGTCCCAGGAAGCACGAAACATGTTGCTTC  
AGAAGCCACAGCTGGCTTATGCTTTGCTGCAAGCACAGGTAGTGTGAGAATTGGGATCCTGAGATTGC  
CCTGAAAATTCTGCACCGCCAGACAAATATCCCAACGCTGATTCAGGCAACCCTCAGCCAGTCCATGTT  
GCCGGACCTGGCTCAGGACCCAATGTTTCAATGAACCAGCAGAATCCTCAGGCCCTCAGGCTCAATCTT  
TGAGTGGAAATGCATGTCATGGCGCACCTCCTATGATGCAAGCTTCTATGCCAGGTGGAGTTCACGACC  
AGTGCAGATGGCAGCTGTGTAGCAGGACCTGGCCCTGGTTTCTTAGCCCTGCAGGAGTCATGCAAGCC  
CAAGTTGGAATGCAAGGAGCTGGACCACTGCCATGGAACGAGGACAAGTCCAAATGCAAGACCCAGAG  
CAGCTATGCAGCGGGGAGCCTTGCCTGCCAACGTCCTCAACTCCTCGTGGTCTTTTAGGTGATGCTCCCA  
TGATCCACGGGGAGGCACTTTAATGACCGTAACGGAGAAGTGAACCTAGAGCTTACCTGGGACCACCA  
CCACCACCTCATCAGGGCCACCGATGCACCACGTTCTGGCCATGAAGGTCGTGGACCACCACCGCATG  
ACATGAGAGGAGGACCGTTAGCTGAGCCAGACCGTTAATGGCAGAGCCAAGAGGACCTATGCTAGATCA  
GAGGGGTCCACCGTTGGATGCCAGAGGTGGAAGAGATCCCCGAGGATTAGATGCACGAGGGATGGAGGCC  
AGAGCCATGGAAGCCAGAGGACTGGAAGCCGTCCTATGGAAGCCGTCCTATGGAAGCCGTCCTATGGA  
AAGCTCGTGCATGGAGGCCGAGCCATGGAGGCTCGTGCATGGAGGCCAGAGGATGGAATACCAGAGG  
CCCAGTACCTGGACCTAGAGGACCTATACCTAGTGAATCCAAGGCCCAATCCAATGAACATGGGGCA  
GTTGTTCCCGAGGATCCAGACAGGTCCAGTGCATGCAGGGAGCAGGCATGCAAGGAGCAAGCATGCAAG  
GTGGAAGCCAGCCTGGTGGCTTTAGTCTGGACAGAGCCAAGTACACCCAGGATCATGAGAAGGCTGC  
TTTGATTATGCAGGTCCTCAACTCACTGCAGACCAGATTGCCATGCTGCCTCCTGAGCAAAGACAGAGC  
ATCCTAATCTAAAGGAACAGATACAGAAATCCACCGGGCACCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001131014

**Insert Size:** 1728 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001131014.1](#), [NP\\_001124486.1](#)

**RefSeq Size:** 4361 bp

**RefSeq ORF:** 1728 bp

**Locus ID:** 683927

**Cytogenetics:** Xq32