

Product datasheet for RN202237

Kazald1 (NM_001033064) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kazald1 (NM_001033064) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Kazald1
Synonyms:	MGC125157
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202237 representing NM_001033064 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCCGAGTGTTCACAGGGCTTCCAGCTAACCATGCTGCACCTTCCTTGGCGCTGCCTTTGGTACTAT
TGTTGCTGATGGTTAGCGCGAGGCCATCCACAGGCCCTGATTACCTGAGGAGAGGCTGGCTGCGGCTGT
AGCCGAGGGCGAGGGTTGTGCTCCCTGCCGGCCAGAAGAGTGCGCCGCGCCGGGGCTGCCTAGCAGGC
CGGTGCGGGATGCGTGGGCTGCTGCTGGGAATGCGCCAACCTGGAGGGCCAGCTCTGCGATCTGGACC
CCAGCGCTAACTTCTACGGGCGCTGCGGCGAGCAGCTCGAATGCAGGCTGGACGCGGGCGGTGACCTGAG
CCGAGGAGAGGTTCCGGAGCCGCTGTGTGTCTGCCGCTCGCAGCGCCCGCTCTGTGGTCCGACGGCCGT
ACCTACGCACAGGTCTGCCGACTGCAGGAGTCCGCCGCTCGCCCGGACGCTAACCTCACTGTGGTGC
ATCCGGGGCCCTGCGAATCGGAGCCCAGATCCTGTGCGAGCCTCACAATATTTGGAACGTGACTGGACA
GGATGTGATCTTTGGCTGTGAGGTGTTTGCCTACCCCATGGCCTCGATTGAGTGGAGGAAAGACGGCTTG
GACATCCAGCTGCCGGGGATGACCCTCATATCTGTGTCAGTTTAGGGTGGACCTCAGAAGTTTGAGG
TAACTGGCTGGCTACAGATTCAAGCTTTGCGCCCTAGTGATGAGGGCACCTATCGCTGCCTTGCCAGGAA
TGCTTTGGCCAGGCTGAGGCCTCTGCTACCCTCACAGTACTCACACCGGACCAGCTGAATGCCACAGGA
TTCTCCAGCTGCAATCCAGGGTTTGAATCCTGAGGAGGAAGCAGAAGATGAAGAGTTGGAAGATTACT
ACTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001033064



[View online »](#)

Insert Size:	915 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033064.1, NP_001028236.1</u>
RefSeq Size:	1645 bp
RefSeq ORF:	915 bp
Locus ID:	293997
Cytogenetics:	1q54