

## Product datasheet for **RN202241**

### Wdr54 (NM\_001109245) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Wdr54 (NM\_001109245) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Wdr54  
**Synonyms:** RGD1561494  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN202241 representing NM\_001109245  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTTCCGCCGAGCGCTCCATCCCACTTCGTGGCTCTGCCGCCCCCTGTCTAACAACTCAGTGTAC  
TGCAGCTTCCAGCCCGTATCTCAGCATTGGGGTGGTTCATGGACCGAGCGCCAGCTTCTCAGCGC  
TGCCCTGAGGGTGTGCCCTGGCCAGCGCCAGCTCCAGTTAAAGAGGGCGCTGGAGTGAGTCCCCCA  
CTATCACTCAGGTCCACTGGTGTGTCTCCCTTCAGAGTACTGCTGGTCTCACCTCACACCGAGGGA  
TACAGATGACGAGTCTGACGGCTCCGTATGGTTTACTGGCATGCCTGGACTCAGGAGATGCTTCTTC  
AGTACAGGCTATGTTGCCGAGGAATTGCTGCCAGCGTCCACTTCATCTGTGTGGGAACGTGGTCTGGT  
CGGATACTGGTCTTTGATATCCCGCTAAGGGTCCCAATATTGTACTCAACGAGGAAGTCTGGGCACC  
AGACACCAATCACAGATATCGCCACTGAGCGTCCAGGGACAGGATGGTGTGCTGACATGGTGACAGC  
AGATGACTCAGGTGTCTGTGTCTGGAGGTGAGGACCTGAGTTCACCTTACTGACCCGAATAGCAGGA  
TTCGGGGTCCCGTGCCCTCTGTGACGCTGTGGCAGGGGATCGTAGCAGCAGGCTATGAAATGGACAAG  
TGCGTCTGTATGACGCTGGCACAGGAGCCCTCCACATCCAGATCAGTGCCACGCCCGGACCATCTGTC  
CCTGGACCTAGTCCAGAGGTGGCAAGCTGTGTCTGACGCTGAGGATACCTTTGTACACATCTGGAAG  
CTGAACAGAAACCCAGAGAGTGGCTCCATTGAGGTGGAACACTGCCATGGCGAATGCATTTCTGACACGC  
AGGTGTGCGGTGCTCGGTTCTGTGACCCAGGGGCGAGCTCCTTTGCTGTGACAGGCTATGACCTTGTGA  
GATCCTGAGATTTGGCAGTGT**TGA**

**ACGGT**ACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001109245



[View online »](#)

<b>Insert Size:</b>	1005 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001109245.1</a></u> , <u><a href="#">NP_001102715.1</a></u>
<b>RefSeq Size:</b>	4354 bp
<b>RefSeq ORF:</b>	1005 bp
<b>Locus ID:</b>	500226
<b>Cytogenetics:</b>	4q34