

Product datasheet for **RN202312**

Vom1r37 (NM_001008923) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vom1r37 (NM_001008923) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Vom1r37
Synonyms:	V1rd6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202312 representing NM_001008923 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGTCTCAAATAAAACCTGAAATTCACAGAGGAAGTGGCTCTTCAGATCCTTTTGCTTTGCCAGT
CTGGGGTTGGGACTGTGGGTAACTGCTCTTTCTTCATAATTTTGTCTAGTTTGGTACTGCTCTCA
ACTGAAGCCCAATCATGTCATTCTATCAAACCTTGTCTGTTGCCAATTTTTATTCTCTTCTCTTGACG
TTTACAAACACTATGGCAGTTTGGGCTCCAAGGAAGCCTCAAAAAACCTCACATGTAACTTACATACT
TCTTTTACATGGTGGCTCGAAGCACAAACATGTGCTCCACCTGTGCCCTGAGCATTTACCAGTTTGTGAC
TCTTGTTCCTGGTAACCTGACTAGGGTCATGCTCAGAGAAATATCTCCAAGGTCATGAAATATTCTTGT
TACAGTTGCTGGTTGTTCAAGTGTCTTAAATAATGCTTACATTCCAATGAATGTCAGTGGTCCACAGAAAA
CAAATAATGACTGATCCTACGGGTAAGTGGGTCTGCTCCATCTCTGGTTTCAGTGTGGCATGAGTTT
CTTGCGGTTTGTCCAGGATACTATTTATCAGCACCATGGTCTGGACCACTGTGTCATGGTCATTAC
TTAAACAGACACCACAGAGAATATGTTACATACACAACCCAGTAAGGACCACAGAGTTATGCTGAGA
CCAGAGCAGCAACACCATCCTGATGCTGGTGGTCACATTTGTGACCTTTTATCTCTAGATTGTATTG
TACTTTTTTTCACATTTCTTTTGTGGACACTCGTTTCTGGTTGAGGCATGTCAAAGAAATTCTGACTGTT
AGCTTCCCCACCATTTCTCCCTTAATATTGATCTTTAGGAATCCTAGGGGTCTTGGCCCTGTACCCCTCG
TTGTGGAAGTCAAAGTCATGTGACTGGGTAGCCAGTCTGGAGCAAGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001008923



[View online »](#)

Insert Size:	963 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008923.1</u> , <u>NP_001008923.1</u>
RefSeq Size:	963 bp
RefSeq ORF:	963 bp
Locus ID:	494272
Cytogenetics:	1q12