

Product datasheet for **RN202322**

Cldn5 (NM_031701) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cldn5 (NM_031701) Rat Untagged Clone
Tag: Tag Free
Symbol: Cldn5
Synonyms: MGC95240
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN202322 representing NM_031701
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGTCTGCAGCGTTGGAAATCTGGGTCTGGTGTGTCTGGTAGGCTGGGTGGGCTTGATCCTGG
CGTGTGGGCTGCCATGTGGCAGGTGACTGCCTTCCTGGACCACAATATCGTGACGGCGCAGACGACTTG
GAAGGGGTGTGGATGTCGTGCGTGGTGCAGAGCACCGGGCACATGCAATGCAAAGTGTATGAGTCTGTG
CTGGCGCTGAGCGCGGAGGTGCAGGCAGCTCGGGCACTACCGTGGCGCTGTGCTGCTGGCGCTTGTGG
CACTCTTTGTTACCTTGACCGCGCTCAGTGCACCACCTGCGTGGCCCCGGGCCAGTTAAGGCACGGGT
GGCACTACGGGAGGAGCGCTTTATGCCCTGTGTGGGCTTCTGGCACTGGTGCCACTCTGCTGGTTCCGCC
AACATCGTAGTCCGGGAGTTCTATGATCCAACGGTGCCGGTGTCTCAGAAGTACGAGCTGGGCGCGGCGC
TGTACATCGGCTGGGCGGCTCCGCACTGCTCATGTGTGGCGGCGGCTCGTGTGCTGTGGCGCTGGGT
TTGCACCGGGCGTCCAGAGTTCAGTTTTCCAGTCAAGTACTCAGCACCAAGGCGAACCACGGCCAACGGC
GATTACGACAAGAAGAACTACGT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1278_h08.zip
Restriction Sites: SgfI-MluI
ACCN: NM_031701
Insert Size: 657 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_031701.2](#), [NP_113889.1](#)

RefSeq Size: 1442 bp

RefSeq ORF: 657 bp

Locus ID: 65131

UniProt ID: [Q9JKD6](#)

Cytogenetics: 11q23

Gene Summary: a component of tight junctions in myelinating Schwann cells [RGD, Feb 2006]