

Product datasheet for **RN202356**

Spata6l (NM_001014165) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spata6l (NM_001014165) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Spata6l
Synonyms:	RGD1310039
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202356 representing NM_001014165 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGCCTCTGGAGGTGGTGGTGGAGCTCCAAATCCGGGCGATTTCTTGCCAGGAGTGTTCTGCCGGACA
AGGAAGTTGTGTACCTTGGAGTCTACCTTCTGAATCAGTACCTGGAGACAGACTGCTTTCCTCTGTGTT
CCCTATTGTGATTCAACAGAGCATGAGGTTTGAAGAGGCAATGGATCCTGGAGCTGTT
GCAGAACTTTAGAAAGCTTCTAACCAGGTTTGAATTAGTACAGCTAGTGTCTCCAGCATGGGAAGAGT
TGGCCTACTATGAAAAAACACCCGAGACTTCTATTCCCTGAGCCCAAGCTGACATCCTCACACCTTGG
CATGCAGAGGGAGGTGCTCATGAAGACGGCTATTTGGTTTCCCGGCAATTGCTCCCAAAGTATGAGTTTCT
ACAAGGACAGCCATCCTAGAGTGTGTGTTCCCGTGTAAGAACAGATTTCTTTGGGAAGAAAGGTGTAAT
TGCAACAGTCTTTCTCAAACTGAACGGACCGGCAAAACAGGAAAAAGGCCAAAGGAGAAAAATTC
TGACCAGTCTCCAAAGGCACACCATTTTGGGGACCTCCCCACAACGACTCCATCTCCACCGGCCACACA
CAGAGGAACCTTGAAAGAGCTTCCGGTTTCTGGGAGAAAGGAAGCCTCCATTTGTGTTAGACATGTGG
ACAGTGGAAACCCCTTTGGTGAAATAATTTGGAGCACCATTACAGAAGTCCAGAAGAAACCTAAGTT
TATAAACTTTGATCTTGCAAAGAAAAGAGCGTTTCTCTTGACAGTCTCGAGGCAACATAAAGGTTGTC
AGAGAGCCGGACGAGCGGATTGTCTTAAGGAGTCAGTCACCGCCACCTCTAAATTCAGGCCATTTTCAGAA
AGCCCTTGCTCAATCACCAAGGGGATGTCGATTTCCACCCGGAAACTACAGTTACCACCTCTCAGCTTTC
CAGCCCTCCCAGCCCTCTGGATCACTCTCTTCTCAAGAAAGGTTCCAGCCTTGTCTGCATCCATGTGG
CAGAAGATCCAGGAGAGGGCATGAGTCTTCTGACACCAACAGAGCTCACCCAAAGGAAGATTCTATCT
CTGAAACAAATTCATCCATGAAAGACCAAGCTACCTTCTGAGGAAACACTTGCTGCGAGAGCACAGATA
TTTT**AA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001014165
Insert Size:	1197 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001014165.1, NP_001014187.1</u>
RefSeq Size:	1534 bp
RefSeq ORF:	1197 bp
Locus ID:	361747
UniProt ID:	<u>Q6AYJ3</u>
Cytogenetics:	1q52