

Product datasheet for **RN202442**

Aldoa (NM_012495) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aldoa (NM_012495) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Aldoa
Synonyms:	Aldo1; RNALDOG5
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN202442 representing NM_012495 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCCCACCCATACCCAGCACTGACCCCGGAGCAGAAGAAGGAGCTGGCTGACATCGCTCACCGAATTG
TAGCTCCGGGCAAGGGCATCCTGGCTGCAGACGAGTCCACTGGAAGCATTGCCAAGCGCTGCAGTCCAT
TGGCACCAGAACACCGAGGAGAACAGGGCCTTCTACCGCAACTGCTGCTGACTGCCGATGACCGTGTG
AATCCCTGCATTGGAGGGGTGATCCTTTCCACGAGACTGTACCAGAAGGCAGATGATGGCCGTCCCT
TCCCCAAGTTATCAAGTCCAAGGGTGGTGTGTGGGCATTAAGGTAGATAAGGGTGTAGTGCCCTGGC
TGAACCAATGGCGAGACCACTACTCAAGGGCTGGACGGGCTGTCTGAGCGCTGTGCCAGTATAAGAAG
GATGGAGCCGACTTTGCCAAGTGGCGCTGTGTGCTAAAGATTGGGGAGCATACTCCCTCGTCCCTCGCCA
TCATGGAAAAATGCCAATGTTCTGGCCCGTTACGCCAGCATCTGCCAGCAGAATGGCATTGTACCCATTGT
GGAGCCTGAAATTCTCCCTGATGGGACCATGACTTGAAGCGCTGCCAGTATGTAAGTGAAGGACTG
GCAGCTGTCTACAAGCTCTGAGTGACCACCATGTCTATCTGGAAGGCACACTGCTGAAGCCCAACATGG
TCACCCCTGGCCATGCTTGACCCAGAAAATTTCCAATGAGGAAATGCCATGGCAACCGTCACAGCATA
TCGTCGAACAGTGCCCCCTGCCGTCCTGGGGTCACTTTCTGTCTGGAGGGCAGAGTGAGGAAGAGGCA
TCCATCAACCTCAATGCTATCAACAAGTGTCCCTGCTGAAGCCATGGGCCTTGACTTTCTCCTATGGCC
GAGCCCTGCAGGCCTCTGCTCTAAAGGCTTGGGGTGGGAAGAAGGAGAACCCTGAAGGCAGCCAGGAGGA
GTACATCAAGCGAGCCCTGGCCAACAGCCTCGCTTGCAAGGAAAGTACACTCCAAGTGGCCAGTCTGGA
GCCGAGCCAGTGAATCTCTCTTATCTCTAACCATGCCTACTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_012495



[View online »](#)

Insert Size:	1095 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012495.2 , NP_036627.1
RefSeq Size:	1465 bp
RefSeq ORF:	1095 bp
Locus ID:	24189
UniProt ID:	P05065
Cytogenetics:	1q36
Gene Summary:	member of the aldolase enzyme family [RGD, Feb 2006] Transcript Variant: This variant (1) represents the longest transcript. All three variants encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.