

Product datasheet for **RN202456**

Pla2g2c (NM_019202) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pla2g2c (NM_019202) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pla2g2c
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN202456 representing NM_019202 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGACCTCTGGTCTCCTCAGGAATGAAGGGCATCGCTGTCTTCTTTATCTTCTGCTGGACAA
CCTCCACCCTCAGCAGCTTCTGGCAGTTCAGAGGATGGTCAAACACATCACGGGGCGCAGCCCTTCTT
CTCCTATTACGGATATGGCTGCTACTGTGGCTTGGGGCCGAGGGATCCCTGTGGACGCCACAGACAGG
TGCTGCTGGGCTCATGACTGTTGCTACCACAAGCTTAAGGAATATGGCTGCCAGCCCATCTGAATGCCT
ATCAGTTTGCCATTGTCAACGGGACCGTGACCTGTGGATGCACCATGGGTGGCGGCTGCTTGTGGGGCA
GAAAGCCTGTGAGTGTGACAACTGTCTGTGTACTGCTTCAAGGAGAACCTGGCCACCTACGAGAAA
TTCAAGCAGCTCTCCCCACCAGCCCCAGTGTGGCAGGGACAACTCCATTGCTAG

ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_019202
Insert Size:	477 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019202.2](#), [NP_062075.1](#)

RefSeq Size: 1757 bp

RefSeq ORF: 477 bp

Locus ID: 29387

UniProt ID: [P39878](#)

Cytogenetics: 5q36

Gene Summary: catalyzes Ca²⁺ dependent hydrolysis of phospholipids including hydrolysis of L-alpha-1-stearoyl-2-arachidonyl phosphatidylinositol [RGD, Feb 2006]