

Product datasheet for **RN202460**

Olr1867 (NM_001101022) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Olr1867 (NM_001101022) Rat Untagged Clone

Tag: Tag Free

Symbol: Olr1867

Synonyms: Gust27; Olr867

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN202460 representing NM_001101022
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGATCCTGAACTGCAACCCATTCTCTGGACTGTTCTGTCCATGTATCTGGTCACAGTGCTTGGGAACT
TGCTCATCATCTGGCTGTCAGCTCTAATTCACATCTCCACAACCTCATGTATTTCTTCTCTCCAATCT
GTCCTTTGTTGACATCTGTTTCATCTCAACCACAATACCAAAAATGCTAGTGAACATACATTCACAGACA
AAAGACATCTCTACATAGAATGCCTTTCACAGGTATATTTTTAACTACTTTTGGTGAATGGATAATT
TTTTACTCACTTTAATGGCCTGTGATCGCTATGTAGCCATCTGCCACCCCTCACTACACTGTAATCAT
GAACCTTCAGCTGTGTGCCCTTCTGATTCTGATGTTTTGGTTAATCATGTTCTGTGTCTCCCTGATTCTAT
GTTCTATTGATGAATGAATTGAATTCTCCAGAGGCACAGAAATCCACATTTCTCTGTGAACTGGCTC
AAGTTCTTAAGGTAGCAATTCTGACACTCATATCAATAATGTCTTCATGTATGTGGTGACTTCCCTACT
AGGACTGATCCCTATGACAGGAATACTTATGTCTTACTCACAGATTGCTTCATCCTTATTAAGATGTCT
TCCTCTGTGAGTAAGTACAAGGCCTTTCCACCTGTGGATCTCACCTCTGTGTGGTCTCTTTATTCTATG
GGTCAGCAACTATAGTTTACTTCTGCTCTTCTGTGCTCCATTCTACACACAAGAAATGATTGCTTCATT
GATGTACACTGTAATCAGCCCCATGCTGAACCCCTTTATCTATAGCCTGAGAAACAAGGATGTAAGGGT
GCCCTTGAAAACTTTTCATCCGAGTTGCCTCTTGCCATTGTGGAGCAAAGACTTTAGACCTAAATTCA
TACTAAACCTGAAAGGCAAAGTTTATA**AA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001101022



[View online »](#)

Insert Size:	939 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001101022.1 , NP_001094492.1
RefSeq Size:	1222 bp
RefSeq ORF:	939 bp
Locus ID:	100036582
UniProt ID:	P34987
Cytogenetics:	8q13
Gene Summary:	Possible olfactory or taste receptor.[UniProtKB/Swiss-Prot Function]