

Product datasheet for **RN202591**

Cox4i1 (NM_017202) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cox4i1 (NM_017202) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Cox4i1
Synonyms:	Cox4; Cox4a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202591 representing NM_017202 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGGCTACCAGGGCACTTAGCCTAATTGGCAAGAGAGCCATTTCTACTTCGGTGTGCCTTCGGGCAC
ATGGGAGTGTGTGAAGAGTGAAGACTATGCTCTCCCATCTTATGTTGATCGGCGTACTACCCCTTGCC
TGATGTGGCCACGTCAAGCTGTCTGCCAGCCAGAAGGCCCTGAAGGAGAAGGAGAAGGCCACTGG
AGCAGCCTTTCCAGGGATGAGAAAGTCCAATTGTACCGCATCCAGTTAACGAGAGCTTCGCTGAGATGA
ACAAGGGCACCAATGAGTGAAGACAGTGGTGGCCTGGCCATGTTCTTCATCGGCTTCACTGCGCTTGT
GCTGATCTGGGAGAAGAGCTACGTGTATGGCCCCATCCCTCATACCTTTGATCGTACTGGTGGCCATG
CAGACCAAGCGGATGCTGGACATGAAGGTCAACCCATTAGGGCTTCTCCGCCAAGTGGGACTACAACA
AGAATGAATGGAAGAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_017202
Insert Size:	510 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017202.1</u> , <u>NP_058898.1</u>
RefSeq Size:	696 bp
RefSeq ORF:	510 bp
Locus ID:	29445
UniProt ID:	<u>P10888</u>
Cytogenetics:	19q12
Gene Summary:	mitochondrial cytochrome C oxidase subunit; enzyme involved in the oxidation of 10-formyltetrahydrofolate [RGD, Feb 2006]