

Product datasheet for **RN202722**

Uts2 (NM_019160) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Uts2 (NM_019160) Rat Untagged Clone
Tag: Tag Free
Symbol: Uts2
Synonyms: Ucn2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN202722 representing NM_019160
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACAGGGTGCCCTTCTGCTGCCTGCTCTTCGTAGGACTCCTGAATCCACTCCTGTCTTTCCCGTCA
CGGACACTGGTGAAATGTCTCTTCAGCTTCCAGTGCTTGAGGAAAATGCTCTTCGGGCTCTGGAGGAGCT
GGAGAGGACTGCCCTCCTGCAGACGCTGCGCCAGACCGTGGGCACAGAAGCAGAGGGAAGCCTAGGCCAG
GCAGATCCAGTGCCGAGACTCCCACTCCAAGGGGAAGCTTGAGGAAGGCTCTCACTGGCAAGATTCTA
ACACTGTACTGAGCCGTCTTTGGCGAGAACCAGGAAACAACGTAAGCAACACGGGACTGCCCCAGAATG
CTTCTGGAAGTACTGCATTTGA****

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_019160

Insert Size: 372 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019160.2](#), [NP_062033.1](#)

RefSeq Size: 527 bp

RefSeq ORF: 372 bp

Locus ID: 29180

UniProt ID: [Q9QZQ4](#)

Cytogenetics: 5q36

Gene Summary: a cyclic neuropeptide; may be involved in the control of neuromuscular functions [RGD, Feb 2006]