

## Product datasheet for **RN202805**

### Cdk1 (NM\_019296) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cdk1 (NM\_019296) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cdk1  
**Synonyms:** Cdc2; Cdc2a  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN202805 representing NM\_019296  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAGGACTATATCAAATAGAGAAAATCGGAGAAGGGACTTATGGTGTGGTATATAAAGGTAGACACA  
 GGACTACTGGCCAGATCGTGGCCATGAAGAAAATCAGACTTGAAGCGAGGAAGAAGGAGTGCCAGTAC  
 GGCAATCCGGGAAATCTCTCTGTAAAAGAACTTCGACATCCAAATATAGTCAGCCTGCAGGATGTACTC  
 ATGCAGGACTCCAGGCTGTATCTCATCTTCGAGTTCCTGTCCATGGATCTCAAGAAGTACCTGGACTCTA  
 TCCTCTCGCCAGTTTATGGATTCTTCGCTCGTTAAGAGTTACTTGTACCAAATCCTCCAGGGGATTGT  
 GTTTTGTCACTCCCGACGAGTTCTTACAGAGACTTGAACCTCAAATCTATTGATTGATGACAAAGGA  
 ACAATCAAACGGCAGATTTTCGGCTTGCAGAGCGTTTGAATACCGATAAGAGTGTACACACACGAGG  
 TAGTGACGCTGTGGTACCGGTCCGAGAGGTGTGCTGGGCTCGCTACTCCACCCCGTTGACAT  
 CTGGAGCATAGGGACCATATTTGCAGAGCTGGCGACCAAGAAGCCGCTTTCCACGGCGACTCAGAGATT  
 GACCAGCTTTCAGGATCTCAGAGCTTGGGCACCCCTAACACGAAGTGTGGCCAGAAGTAGAGTCCC  
 TGCAGGACTACAAGAACACTTTCCCAAGTGAAGCCGGGAGCCTCGGTCCCACGTCAAGAACCTGGA  
 TGAACACGGCTTGGACTTGTCTCGAAAATGCTGGTCTATGATCCAGCCAAACGAATCTCTGGAAAATG  
 GCCCTAAGCACCCATACTTTGACGACTTGGACAATCAGATCAAGAAGATG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_019296  
**Insert Size:** 894 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_019296.1</a> , <a href="#">NP_062169.1</a>
<b>RefSeq Size:</b>	1184 bp
<b>RefSeq ORF:</b>	894 bp
<b>Locus ID:</b>	54237
<b>UniProt ID:</b>	<a href="#">P39951</a>
<b>Cytogenetics:</b>	20p11
<b>Gene Summary:</b>	putative cell cycle control protein; homologous to <i>S. pombe</i> Cdc2 [RGD, Feb 2006]