

## Product datasheet for **RN202854**

### Tgif2 (NM\_001134983) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tgif2 (NM_001134983) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tgif2
Synonyms:	RGD1560115; RGD1564927; Tgif2l1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202854 representing NM_001134983 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCGGACAGCGATCTAGGCGAGGATGAAGTCTTCTGTCACTGACTGGCAAGAGGAAGCGCAGGGGGA  
ACCTGCCCAAGGAGTCAGTAAAGATTCTCCGAGACTGGCTGTATCTGCACCGCTACAACGCCTACCCCTC  
TGAGCAGGAGAAGCTAAGCCTCTCTGGACAGACCAACCTCTCCGTGCTGCAGATATGTAAGTGGTTCATC  
AATGCCCGCCCGCCTTCTCCCTGACATGCTTCGGAAGGATGGCAAAGACCCTAATCAGTTCACCATCT  
CCCGGCGGGGGTAAGGCCTCAGATGTGGCTCTCCCCAGGGCAGCAGCCCTCTCTGTGGCTGTGTC  
TGTTCCAGCCCCACCAACATGCTCTCCTGTCTGTGTTCATGCCACTCCACTCAGGCCAAGGTGAA  
AAGCCAGCAGCGTCTTCCCACAAGTGGAGCTGGAGTCCCCCAAGGCCTGGTGACCCCTGGGAGCACAC  
TCACCCTGCTGACTAGGGCAGAGGCTGGAAGCCCCACAGGTGGACTCTTCAATACGCCACCCACACC  
CCCAGAGCAGGACAAGGACGACTTCAGCAGCTTCAGTTGCTGGTGGAGGTGGCGCTGCAGAGGGCTGCT  
GAGATGGAGCTTCAGAAGCAGCAAGATCCAGCGCCACCCTTACTACACTCCGCTCCCTCTAGTCTCAG  
AAAATGCCAAGTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001134983
Insert Size:	714 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001134983.1</a></u> , <u><a href="#">NP_001128455.1</a></u>
<b>RefSeq Size:</b>	2908 bp
<b>RefSeq ORF:</b>	714 bp
<b>Locus ID:</b>	499929
<b>Cytogenetics:</b>	3q42