

## Product datasheet for **RN202897**

### Chst14 (NM\_001109639) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Chst14 (NM\_001109639) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Chst14  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN202897 representing NM\_001109639  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGC**C

ATGTTTCCCCGCCCTCTGACCCCACTAGCTGCTCCGAAAAGCGCCAAGACCCTGGCCCGCAGCCGAAGC  
GGCCCCGTGGGCAGGGCCCGGGCTGGGATCGGGGGCCGCCCTGCTACTGCCGTCCATGCTGATGTT  
CGCAGTAATCGTGGCGTCCAGCGGACTGCTCATGATCGAGCGGGGCATCCTATCGGAGATGAAACCC  
CTTCCCCTGCACCCTCCAGCCACAAAGGCGCGCCTGGAGCGGGACAGATCCTAAGCCTAGAGGCCTAT  
CTTTCGATGCTGAGGACTCGGACTTGAAGTGCGGGAGGACATCCGAAACCGGACCTTGAGGGCCGTGTG  
CGGACAACCAGGCATGCCCGGGACCCCTGGGACTTGCCGGTGGGACAGCGGCGCACCCCTACTACGCCAC  
ATCCTCGTAAGTGACCGTACCGCTTCTCTACTGCTACGTCCCCAAAGTGGCCTGCTCTAACTGGAAC  
GCGTGTGAAGGTGCTGGCAGGCGTCTGAACAACGTGGATGTCCGCCTCAAGATGGACCACCGCAGTGA  
CTTGGTGTCTGGCAGACCTGCGGCCTGAGGAGATTCGCTACCGTCTGCAGCACTACTCAAGTTCTGT  
TTTGTGCGAGACCCCTTGAACGCCTCCTGTCTGCTTACCGCAACAAGTTGGAGAGATCCGAGAGTACC  
AGCAGCGCTATGGGGCCGAAATTGTGAGGCGCTACAGGGCTGGAGCTGGTCCCAGCCCTGCGGGGACGA  
CGTTACCTTCCCAGAGTTCCTGAGATATCTGGTGGATGAAGATCCTGAGCATATGAATGAGCATTGGATG  
CCTGTGTACCACCTGTGCCAACCATGTGCTGTGCACTATGACTTCGTGGGTTCCATGAGAGGCTGGAAG  
CTGATGCCAACAGGTGCTGGAGTGGGTGCGGGCCCCATCCCATGTCCGGTCCCAGCTCGCCAGGCCTG  
GTACCGTCCAGCCAGCCAGAAAGTCTACATTACCACCTGTGCAATGCCCGCGGGCCCTGCTTCAAGAT  
GTGCTACCTAAGTATATCCTAGACTTCTCCCTTTGCTTACCCACTGCCAATGTCACCAAGGAAGCCT  
GTCGC**TAG**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-RsrII  
**ACCN:** NM\_001109639



[View online »](#)

<b>Insert Size:</b>	1128 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001109639.1, NP_001103109.1</u>
<b>RefSeq Size:</b>	2058 bp
<b>RefSeq ORF:</b>	1128 bp
<b>Locus ID:</b>	691394
<b>Cytogenetics:</b>	3q35