

Product datasheet for **RN202980**

Csn2 (NM_017120) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Csn2 (NM_017120) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Csn2
Synonyms:	Csnb
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN202980 representing NM_017120 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGTCTTCATCCTTGCCTGCCTTGTGGCACTTGCTCTTGCAAGGGAGAAGGATGCATTCACCTGTGT
CCTCTGAGACTGGTAGTATTTCCAGTGAGGAATCCGTTGAACATATCAATGAGAACTTCAGAAGGTGAA
ACTCATGGGACAGGTGCAGTCCGAGGATGTGCTCCAGAATAAATTTCACTCCGGCATTCAAGTCAACCC
CAGGCCATTCCATATGCCAGACCATCTCTGCAGTCCCATTCCACAAAACATCCAGCCTATTGCTCAAC
CCCCTGTGGTGCCAACTGTTGGGCTATCATTTCTCCTGAACTGGAATCCTTCTTAAAGCTAAAGCCAC
TGTCTTCCCAAGCACAACAGATGCCCTTCTTAACTCTGAAACCGTGCCTCCGCCTTTAACTCTCAA
ATCCCCAGTCTTGATCTTGCTAATCTGCACCTTCTCAGTCTCCAGCCAGCTCCAGGCACAAAATTGTG
AGGCTTTTCTCAGACTCCCGGGTGGTTTCTTCTCAGCCCCAGCTGTCTTCTCAGTCCAAAAGTCA
GTACCTTGTGCAGCAACTAGCACCCCTTCCAACAAGGTATGCCTGTTCAAGACCTTCTGCAGTACCTA
GATCTTCTGCTTAACCCACCCTCCAGTTCCTTGCCACTCAACAACCTTATTCCACTTCTGTCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_017120
Insert Size:	696 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017120.3</u> , <u>NP_058816.2</u>
RefSeq Size:	1130 bp
RefSeq ORF:	696 bp
Locus ID:	29173
UniProt ID:	<u>P02665</u>
Cytogenetics:	14p21
Gene Summary:	expression is induced by pseudophosphorylated prolactin in mammary gland [RGD, Feb 2006]