

## Product datasheet for **RN203075**

### Socs4 (NM\_001107256) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Socs4 (NM\_001107256) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Socs4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN203075 representing NM\_001107256  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTGAAAACAATAGTAAAAATGTAGATGTGCGGCCTAAAACAAGTCGGAGTAGAAGTGCTGACAGGA  
AGGACGGTTATGTGTGGAGTGGAAAGAAGTTATCTTGGTCCAGAAAGAGTGAGAGCTGTTCTGAGTCTGA  
GGCCGTAGGTACTCTTGAGAGCGTTGAAATCCCTCTGAGGAGCCAGGAAAGGGAACCTAGCTGTTCTGTC  
ATTGAGTTGGACTTGGATCATTCCCTGTGGGCATAGATTTTTAGGCCGATCCCTTAAACAGAAACTGCAAG  
ATGCCGTGGGGCAGTGTTTTCTATAAGGAACTGTGGTCTCGCCACTCTCCAGGGCTTCCATCTAAAAG  
AACGATTCACATCAGTGAGCTCATGCTAGACAAGTGCCCTTTCCACCTCCATCAGACTTAGCCTTAGG  
TGGCATTATTAACCGGCACACTGTTGCTCTGAGTTCAGTTCAGATGAATGGTGAGCACAGACCTGT  
GTGAGAGGAGGCCGAGAGATGCTCAGCTGAAACGAAGAAGCACAGAAGACGACATACCCTGTTTCTACA  
TACCAGTGTCAGCCGTGTGCATAACTACCAACAGTGCTTCATGTGCAGGTGGTCATGTGACCGGTTCT  
CTGATGAACTTGGTCACCGACAATAGCATGGAAGACGGTGATATGGACTCAGAGGATGAAATTATAACAC  
TGTGCACAAGCTCCAGGAAAAGGAACAAGCCAGATGGGAAACGGATGACGGCATCCTGCAATTGGAGGC  
ACCTCCTAAGTTCCACACCCAGATTGATTACGTGCACTGCCTTGCCAGACCTCCTCCAGATCAGTAAC  
AACCCGTGCTACTGGGGGTTATGGACAAATACGCAGCCGAGGCTCTGCTAGAAGGAAAGCCAGAGGGCA  
CCTTTTTACTTCGAGATTCAGCCAGGAAGATTATTTATCTCTGTTAGTTTTAGGCCCTACAGTCGTTT  
TCTTCACGCTAGAATTGAACAGTGGAATCATAACTTCAGCTTTGACGCCACGATCCCTGTGCTTCCAT  
TCTCCGATATTACTGGGCTCCTAGAGCATTATAAGGACCAAGTGCCCTGCATGTTCTTTGAGCCACTTT  
TATCCACTCCGTTAATCCGGACTTTCCCTTTTCTTGCAGCACATTTGCAGGACAGTATTTGTAACTG  
TACAACGTACGATGGCATCGATGCCCTTCAATTCCTTCTCTATGAAATTGTACCTGAAGGAATACCAT  
TACAAATCAAAGTTAGGCTACTCAGGATTGACGTGCCAGAGCCACAGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_001107256
<b>Insert Size:</b>	1311 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001107256.2</a></u> , <u><a href="#">NP_001100726.1</a></u>
<b>RefSeq Size:</b>	2037 bp
<b>RefSeq ORF:</b>	1311 bp
<b>Locus ID:</b>	305828
<b>Cytogenetics:</b>	15p14