

## Product datasheet for **RN203187**

### Cryge (NM\_173289) Rat Untagged Clone

#### Product data:

|                           |                                                                                 |
|---------------------------|---------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                             |
| Product Name:             | Cryge (NM_173289) Rat Untagged Clone                                            |
| Tag:                      | Tag Free                                                                        |
| Symbol:                   | Cryge                                                                           |
| Synonyms:                 | Cryg5; Len                                                                      |
| Mammalian Cell Selection: | Neomycin                                                                        |
| Vector:                   | pCMV6-Entry (PS100001)                                                          |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                            |
| Fully Sequenced ORF:      | >RN203187 representing NM_173289<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGGAAGATCACCTTCTATGAGGACCGCGCTTCCAGGGCCGCACTATGAGTGCAGCACAGACCACT  
CCAACCTGCAGCCCTACTTCAGCCGCTGCAACTCTGTGCGCTGGACAGTGGCTGCTGGATGCTCTATGA  
GCAGCCCAACTTCACAGGCTGCCAGTATTTCTTCGTCGCGGGGACTACCCTGACTACCAGCAGTGGATG  
GGTTTCAGCGACTCTGTCCGCTCCTGCCGCTCATCCCCACTCCAGCTCTACAGAATCAGGATCTACG  
AGCGAGAGGACTACAGAGGCCAGATGGTGGAGATCACAGACGACTGCCCCACCTGCAGGACCGCTCCA  
CTTCAGTGACTTCCACTCTTCCACGTGATGGAGGGCTACTGGGTCTCTATGAGATGCCAACTACCGG  
GGGCGGCAGTACCTGCTGAGGCCTGGAGAATACAGGCGCTACCACGACTGGGCGCCATGAATGCCAGG  
TAGGCTCTCTGAGGAGAATCATGGATTTCTAT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

|                    |                                                                                                                                                                                                                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Restriction Sites: | SgfI-MluI                                                                                                                                                                                                                                                                                                                                      |
| ACCN:              | NM_173289                                                                                                                                                                                                                                                                                                                                      |
| Insert Size:       | 525 bp                                                                                                                                                                                                                                                                                                                                         |
| OTI Disclaimer:    | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |


[View online »](#)

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Components:            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| Reconstitution Method: | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| RefSeq:                | <u>NM_173289.1, NP_775411.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| RefSeq Size:           | 633 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:            | 525 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:              | 24279                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| UniProt ID:            | <u>P02528</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Cytogenetics:          | 9q32                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene Summary:          | putative protein of the lens of the eye [RGD, Feb 2006]                                                                                                                                                                                                                                                                                                                                                                                                                                                           |