

Product datasheet for RN203192

Pdpdf (NM_001009316) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pdpdf (NM_001009316) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pdpdf
Synonyms:	RGD1310660
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN203192 representing NM_001009316 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGCTATCCCTCCAGCGGCTCGCTCGTGGCTACCCATGACTACTATCGGCGACGTCTGGGCTCCT
CGTCCAGCAACAGCTCCGGCGGAAGTGCAGAGTACCCGGGGACGCCGTGCCCACTCCCGGGTCTTCC
CAAAGCTGACTCGGGCCACTGGTGGCCAGCTTTTTCGGGAAATCCACCCTCCATTATGGCCGA
GTGTTGGAGTCTCTGAGCACTCGGCAGAATCCCCCAGGCCTCCAGAAGCCGATCTCCTGTGGCCTGG
CTCCTGAAACCATGAAGAAGCAGCCAGTCATGCATCCTAGCCAGACCAACAGCAGGGCCCGTCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001009316
Insert Size:	348 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001009316.1](#), [NP_001009316.1](#)

RefSeq Size: 685 bp

RefSeq ORF: 348 bp

Locus ID: 296470

UniProt ID: [Q5PR01](#)

Cytogenetics: 3q43

Gene Summary: Probable regulator of exocrine pancreas development.[UniProtKB/Swiss-Prot Function]